


Mepitel[®] One

THE WOUND CONTACT LAYER WITH ONE SIDED SAFETAC[®] TECHNOLOGY

- Minimises pain to the patient and trauma to the wound at dressing changes^{1,2}
- Can remain in place for up to 14 days³⁻⁷ which allows cost-effective^{1,8} and undisturbed wound healing




Polyurethane net

- Open mesh structure allows free passage of exudate into secondary dressing
- High transparency allows improved wound visibility
- Allows for application of topical treatment


Safetac[®] layer

- Reduces pain and trauma before, during and after dressing changes⁹
- Does not adhere to the moist wound bed but to dry skin only
- Seals the wound margins and reduces the risk of maceration¹⁰

HOW TO USE MEPITEL® ONE



1. Clean the wound area. Remove the release film.



2. Apply Mepitel® One to the wound.



3. Apply outer absorbent dressing, e.g. Mesorb.

HOW MEPITEL® ONE WORKS

Mepitel One can be left in place for up to 14 days³⁻⁷, depending on the condition of the wound which reduces the necessity for frequent primary dressing changes. The porous structure of Mepitel One allows exudate to pass into a secondary dressing. The Safetac layer prevents the outer dressing from sticking to the wound and ensures atraumatic dressing changes. The Safetac layer seals around the wound edges, preventing the exudate from leaking onto the surrounding skin, thus minimising the risk of maceration.

BENEFITS OF MEPITEL® ONE

- Minimises pain and trauma at dressing changes.
- Can remain in place for up to 14 days which in turn ensures undisturbed wound healing
- Enables less frequent dressing changes
- Minimises the risk of maceration
- Maintained product properties over time – leaves no residues and does not dry out
- Transparent for easy wound inspection during application and during wear
- Conforms well to body contours, promoting patient comfort during wear
- Stays in place
- Low potential for skin irritation and allergy



References:

1. Gotschall CS, et al. Prospective, randomized study of the efficacy of Mepitel on children with partial-thickness scalds. *Journal of Burn Care & Rehabilitation* 1998;19(4):279-283.
2. Bugmann Ph, et al. A silicone-coated nylon dressing reduces healing time in burned paediatric patients in comparison with standard sulfadiazine treatment: a prospective randomized trial. *Burns* 1998;24:609-612.
3. Eagle M. Use of non-adherent silicone dressing Mepitel to meet client centered needs in chronic non-healing wounds. *Clinical Report Mölnlycke Health Care* (1998).
4. Taylor R. Use of a silicone net dressing in severe mycosis fungoides. *JoWC*, Vol 8, No 9 (1999), p.429-430.
5. Young T. Fungating wounds: their diagnosis and management. *Community nurse*, 5, No 10 (1999), p. 53-54.
6. Marconi R, Laverda F, Trevisan G. Poster Presentation European Wound Management Association, Helsinki, Finland, 2009.
7. Barraziol R, Fracalenta E, Schiavon M. Poster Presentation European Wound Management Association, Helsinki, Finland, 2009.
8. Rippon M, Davies P, White R, Bosanquet N. Cost implications of using an atraumatic dressing in the treatment of acute wounds. *JoWC*, vol 17, No 5 (2008), p. 224-7.
9. Dykes PJ et al. Effects of adhesive dressings on the stratum corneum of the skin. *J Wound Care* 2001; 10(2):7-10.
10. Dykes PJ. The effect of adhesive dressing edges on cutaneous irritancy and skin barrier function. *J Wound Care* 2007, 16(3):97-100.

PRECAUTIONS:

- The wound should be inspected for signs of infection according to clinical practice. Consult a healthcare professional for the appropriate medical treatment.
- If Mepitel One is used on Epidermolysis Bullosa patients, employ extra surveillance at dressing changes. Adhesion level of Mepitel One is somewhat higher than of Mepitel.
- When used on bleeding wounds or wounds with high viscosity exudate, Mepitel One should be covered with a moist absorbent dressing pad.
- When Mepitel One is used for the fixation of skin grafts, the dressing should not be changed before the fifth day post application.

AREAS OF USE

Mepitel One is designed for a wide range of wounds such as painful wounds, skin tears, skin abrasions, surgical incisions, partial thickness burns, traumatic wounds, blistering, lacerations, partial and full thickness grafts, radiation skin, leg and foot ulcers. It can also be used as a protective layer on non-exuding wounds and on areas with fragile skin.

WARNINGS

When Mepitel One is used on burns treated with meshed grafts or after facial resurfacing imprints can occur if the product is not used properly.

MEPITEL ONE ORDERING INFORMATION

Art. No	Size cm	Pieces Per Inner	PIP Code	NHS Code
289170	6 x 7	5	350-6888	EKH037
289270	9 x 10	5	350-6896	EKH038
289470	13 x 15	5	350-6904	EKH039
289670	24 x 27.5	5	350-6912	EKH040



Mölnlycke Health Care Ltd, Arenson Centre, Arenson Way, Dunstable, Bedfordshire LU5 5UL, UK. Tel. +44 (0) 870 6060766, www.molnlycke.co.uk

The Mölnlycke Health Care name and logo, Safetac® and Mepitel® One are registered trademarks of Mölnlycke Health Care AB. ©2009 Mölnlycke Health Care. All rights reserved.