

BARRIER® EasyWarm® active self-warming blanket

- Simple and quick to set-up and remove
- Portable across the patient journey due to no extra devices being needed
- Easily available to more patients thanks to no need for additional equipment
- Noiseless
- No air turbulence thus avoiding spread of particles and microorganisms



BARRIER will help you treat more patients to better care

All BARRIER products from Mölnlycke Health Care are designed to offer the optimal combination of protection, functionality and comfort. Our innovative and procedure specific solutions come complete with logistic advantages and expert training and support. These elements together will create a ripple effect of savings in time, effort and costs - helping you treat more patients to better care.

User instructions



Tear package open using pre-cut tab located in top right corner. Remove the blanket from the outer protective package to expose BARRIER® EasyWarm® to air.



Open and unfold BARRIER EasyWarm completely and do not fold it over itself. Blanket will self-warm to an average of 40°C in 30 min once exposed to air, and keep its temperature up to 12 hours.



Place BARRIER EasyWarm with pads up. No extra weight should be added on top of warming pads.



During the pre- and post operative period and during surgical interventions on hand, arms and head place BARRIER EasyWarm over the entire body.



For surgical intervention in the lower abdominal area and lower extremities place BARRIER EasyWarm over the upper extremities and chest, with the cutout of the blanket at head level.



For surgical interventions above the genital area place BARRIER EasyWarm over the lower part of the body.

Instructions for use

- To be handled under supervision of healthcare provider.
- Avoid using in an oxygen rich environment.
- BARRIER EasyWarm should be removed from patient if an x-Ray or MRI is needed since the warming pads are not translucent.
- It is recommended to avoid direct contact with bruising, swelling, frostbite or desensitized skin.
- After use, dispose BARRIER EasyWarm with regular waste. Possible to incinerate.
- Warming pads should be placed upwards*

*



Ref. No.	Description	Trp
629995-00	<p>BARRIER EasyWarm</p> <p>Product: 12 air activated warming pads (13cmX10cm) integrated into pockets of the top layer. Blanket size is 152cmX92cm.</p> <p>Material: Several layers of SMMMS (polypropylene). The layers and pockets are sewn together with Polyester/cotton threads. The warming pads are made up of air permeable bags that contain active coal, clay, salt and iron powder.</p> <p>Packaging: Vacuum packed into an airtight polymer bag. The bags (10pcs) are further packed into a polyethylene bag lined transport box with bubble wrap.</p>	10

Precautions

BARRIER EasyWarm is not designed to be used on infants. The blanket should not be folded over itself. If prolonged hypotension can be expected, use of the blanket should be avoided. The blanket should not be used in cases of severe peripheral vascular disease.

