




Transfusion Sets for Blood or Blood Components


- DEHP-free Single and Double Chamber Sets with 200µm filtration
- 3mm and 4mm ID tubing
- 20 drops per minute
- White flotation ball in Double Chamber Sets
- Available with 3-way connector and extension tubing
- Trauma pump variant to assist rapid emergency infusions
- Dual spike option

Single Chamber Blood Sets (3mm ID)



7009DF3MM


Single chamber Blood Set with 200µm filter, 190cm tubing (3mm ID), rotating luer connector



7009DF3MMXMY


Single Chamber Blood Set with 200µm filter, 190cm tubing, (3mm ID), 3-way connector, needle free access port, 120cm extension tubing, rotating luer connector

Single Chamber Blood Sets (4mm ID)




7009DF4SML

Single Chamber Blood Set with 200µm filter, 180cm tubing (4mm ID), rotating luer connector



7009DF4XSML


Single Chamber Blood Set with 200µm filter, 180cm tubing (4mm ID), 3-way connector, 120cm extension tubing, rotating luer connector



7009DF4XMY SML


Single Chamber Blood Set with 200µm filter, 180cm tubing (4mm ID), 3-way connector, needle free access port, 120cm extension tubing, rotating luer connector

Double chamber Blood Sets (4mm ID)




7017DF4

Double Chamber Blood Set with 200µm filter, 142cm tubing (4mm ID), rotating luer connector




7017DF4X

Double Chamber Blood Set with 200µm filter, 142cm tubing (4mm ID), 3-way connector, 142cm extension tubing, rotating luer connector




7017DF4XMY

Double Chamber Blood Set with 200µm filter, 142cm tubing (4mm ID), 3-way connector, needle free access port, 142cm extension tubing, rotating luer connector



7017DF4D

Dual Spiked, Double Chamber Blood Set with 200µm filter, 142cm tubing (4mm ID), rotating luer connector



7017DF4TP

Double Chamber Blood Set with 200µm filter, trauma pump, 142cm tubing (4mm ID), rotating luer connector



Single Chamber Blood Sets with 3mm or 4mm ID tubing



Available with 3-way connector and extension tubing



Needle free access port as optional extra tubing

| Product code | Product description | Units per box | NHS Supply Chain code |
|---|---|---------------|-----------------------|
| Single Chamber Blood Sets (3mm ID) | | | |
| 7009DF3MM | Single Chamber Blood Set with 200µm filter, 190cm tubing (3mm ID), rotating luer connector | 50 | – |
| 7009DF3MMXMY | Single Chamber Blood Set with 200µm filter, 190cm tubing, (3mm ID), 3-way connector, needle free access port, 120cm extension tubing, rotating luer connector | 50 | – |
| Single chamber Blood Sets (4mm ID) | | | |
| 7009DF4SML | Single Chamber Blood Set with 200µm filter, 180cm tubing (4mm ID), rotating luer connector | 50 | – |
| 7009DF4XSML | Single Chamber Blood Set with 200µm filter, 180cm tubing (4mm ID), 3-way connector, 120cm extension tubing, rotating luer connector | 50 | – |
| 7009DF4XMY SML | Single Chamber Blood Set with 200µm filter, 180cm tubing (4mm ID), 3-way connector, needle free access port, 120cm extension tubing, rotating luer connector | 50 | – |
| Double chamber Blood Sets (4mm ID) | | | |
| 7017DF4 | Double Chamber Blood Set with 200µm filter, 142cm tubing (4mm ID), rotating luer connector | 20 | FSB1880 |
| 7017DF4X | Double Chamber Blood Set with 200µm filter, 142cm tubing (4mm ID), 3-way connector, 142cm extension tubing, rotating luer connector | 20 | – |
| 7017DF4XMY | Double Chamber Blood Set with 200µm filter, 142cm tubing (4mm ID), 3-way connector, needle free access port, 142cm extension tubing, rotating luer connector | 20 | – |
| 7009DF4D | Dual spiked, Double Chamber Blood Set with 200µm filter, 142cm tubing (4mm ID), rotating luer connector | 20 | – |
| 7009DF4TP | Double Chamber Blood Set with 200µm filter, trauma pump, 142cm tubing (4mm ID), rotating luer connector | 20 | – |