



Advanced Leg Ulcer Treatment

The Altipress® 40 Leg Ulcer Hosiery Kit offers a practical and discreet alternative to bulky bandage systems. The kit consists of two stockings that are worn one on top of the other. The liner offers 10mmHg and the overstocking is a British Standard Class III stocking. Together, they provide the 40mmHg pressure at the ankle recommended for leg ulcer healing.

Altipress® Liners are available separately in packs and provide 10mmHg at the ankle. They are used for palliative care, maintenance therapy and for those intolerant of higher levels of compression. They can also be layered to provide higher levels of compression.

Benefits:

- Five stock sizes in 3 lengths
- Made to Measure available
- Latex Free
- Cost effective
- Easy to apply
- Comfortable fit
- Ideal for self-caring patients
- Made in the UK

Indications:

- Treatment of and management of venous leg ulcers
- Prevention of VLU recurrence prevention
- Patients not suitable for or tolerant of compression bandages

Contra-indications:

- Peripheral Vascular Disease (PVD). This may be confirmed with an ABPI (Doppler) reading of less than 0.8
- Acute or unstable Heart Failure
- Ischaemic Heart Disease (IHD), also called Coronary Artery Disease
- Untreated Cellulitis or Phlebitis of the Lower Limb

Proceed with caution in those patients who have:

- Diabetes
- Rheumatoid Arthritis
- Heavily exuding or large wounds
- Irregular shaped limbs

These patients may require referral to a specialist for advice, or they may require an additional assessment and/or close supervision under the care of an appropriately trained healthcare professional in line with local policies/pathways





Stock Range

Kits:

- Contain 2 x liners & 1 x outer stocking
- Liners available in White
- Outer stockings available in Soft Beige

Liner Packs:

- Contain 3 x liners
- Available in White or Soft Beige


FP10 / GP10 Ordering Codes (PIP Codes)

Kit	Altipress 40* Hosiery Kit	Altipress 40* Hosiery Kit	Altipress 40* Hosiery Kit	Altipress* Liner Pack					
Length	Short	Regular	Long	Short	Short	Regular	Regular	Long	Long
Colour	Soft Beige	Soft Beige	Soft Beige	White	Soft Beige	White	Soft Beige	White	Soft Beige
Small	337-1325	337-1374	337-1424	337-1473	341-8506	337-1523	341-8662	337-1572	341-8598
Medium	337-1333	337-1382	337-1432	337-1481	341-8696	337-1531	341-8654	337-1580	341-8563
Large	337-1341	337-1390	337-1440	337-1499	341-8688	337-1549	341-8647	337-1598	341-8530
X Large	337-1358	337-1408	337-1457	337-1507	341-8670	337-1556	341-8621	337-1606	341-8522
XX Large	337-1366	337-1416	337-1465	337-1515	341-8613	337-1564	341-8605	337-1614	341-5414

Measuring and Fitting Guide

	B: Calf	C: Ankle
Small	32-38	19-25.5
Medium	34-40	21-27.5
Large	36-42	23-29.5
X Large	38-45	25-31.5
XX Large	40.5-50	26.5-34

	Leg Length
Short	34-38
Regular	38-42
Long	42-46



Measurements in cms

Made to Measure Range

Kits:

- Contain 2 x liners & 1 x outer stocking
- Liners available in White
- Outer stockings available in Black or Soft Beige

Liner Packs:

- Contains 2 x liners
- Available in White



Made to Measure Order forms are available from
www.altimed.co.uk or by contacting Customer Services



01509 501 720



01509 501 721



enquiries@altimed.co.uk