

# AMSafe®

## Pre-Filled Normal Saline Flush Syringe

0.9% Sodium Chloride

### Features

- Provides all convenient size syringes and volumes
- DEHP-free, PVC-free and latex-free
- Sterile, non-pyrogenic and non-toxic
- Standard closure cap and sterile field configurations
- Reduced rebound
- Eliminates needles and vials
- Eliminates potential for cross-contamination
- Complies with all PICC and INS standards
- Compatible with ISO standard needle-free valves
- Properly licensed, insured, fully FDA cleared and CE marked



### Sterile Field Ready

Product Code	Product description	Units per box	NHSSC code
IVF1203SS	AMSafe® 3mL Pre-filled 0.9% Sodium Chloride Flush Syringe (Sterile Field Ready)	30	FWC588
IVF1205SS	AMSafe® 5mL Pre-filled 0.9% Sodium Chloride Flush Syringe (Sterile Field Ready)	30	FWC589
IVF1210SS	AMSafe® 10mL Pre-filled 0.9% Sodium Chloride Flush Syringe (Sterile Field Ready)	30	FWC590
IVF2020SS	AMSafe® 20mL Pre-filled 0.9% Sodium Chloride Flush Syringe (Sterile Field Ready)	20	FWC591

### Sterile Fluid Pathway

Product Code	Product description	Units per box	NHSSC code
IVF1203	AMSafe® 3mL Pre-filled 0.9% Sodium Chloride Flush Syringe (Sterile Fluid and Path)	30	FWC665
IVF1205	AMSafe® 5mL Pre-filled 0.9% Sodium Chloride Flush Syringe (Sterile Fluid and Path)	30	FWC664
IVF1210	AMSafe® 10mL Pre-filled 0.9% Sodium Chloride Flush Syringe (Sterile Fluid and Path)	30	FWC663
IVF2020	AMSafe® 20mL Pre-filled 0.9% Sodium Chloride Flush Syringe (Sterile Fluid and Path)	20	FWC662