

Neutral Displacement

Needle-free connector



**A1010N**

NeutrArt® Needle-free Valve Neutral

Priming Volume	Length	Flow Rate
0.04 ml	2.5 cm	160 ml/min

AU1010N

NeutrArt® Needle-free Valve Neutral with Extension Line

Priming Volume	Length	Flow Rate
0.19 ml	10 cm	70 ml/min

AUWB1010N

NeutrArt® Needle-free Valve With 4.0 mm OD Extension Line

Priming Volume	Length	Flow Rate
0.46 ml	10 cm	90 ml/min

A2010N

NeutrArt® Double Lumen Needlefree Valve Neutral

Priming Volume	Length	Flow Rate
0.42 ml	10 cm	70 ml/min

A3010N

NeutrArt® Triple Lumen Needle-free Valve Neutral

Priming Volume	Length	Flow Rate
0.57 ml	10 cm	70 ml/min

AU1010NWOR

NeutrArt® Needle-free Valve Neutral with Extension Line, Luer Slip

Priming Volume	Length	Flow Rate
0.19 ml	10 cm	70 ml/min

A2010NWOR

NeutrArt® Double Lumen Needle-free Valve Neutral, Luer Slip

Priming Volume	Length	Flow Rate
0.42 ml	10 cm	70 ml/min

A3010NWOR

NeutrArt® Triple Lumen Needle-free Valve Neutral, Luer Slip

Priming Volume	Length	Flow Rate
0.57 ml	10 cm	70 ml/min

ATU3010N

NeutrArt® Needle-free Valve with Three Way Stopcock, Neutral

Priming Volume	Length	Flow Rate
0.48 ml	16.5 cm	60 ml/min

**A1010N-PH**

NeutrArt® Needle-free Valve Neutral, Photoprotective

Priming Volume	Length	Flow Rate
0.04 ml	2.5 cm	160 ml/min

AU1010N-PH

NeutrArt® Needle-free Valve Neutral with Extension Line, Photoprotective

Priming Volume	Length	Flow Rate
0.19 ml	10 cm	70 ml/min

AUWB1010N-PH

NeutrArt® Needle-free Valve With 4.0 mm OD Extension Line, Photoprotective

Priming Volume	Length	Flow Rate
0.46 ml	10 cm	91 ml/min

A2010N-PH

NeutrArt® Double Lumen Needlefree Valve Neutral, Photoprotective

Priming Volume	Length	Flow Rate
0.42 ml	10 cm	70 ml/min

A3010N-PH

NeutrArt® Triple Lumen Needle-free Valve Neutral, Photoprotective

Priming Volume	Length	Flow Rate
0.57 ml	10 cm	70 ml/min

AU1010NWOR-PH

NeutrArt® Needle-free Valve Neutral with Extension Line, Luer Slip, Photoprotective

Priming Volume	Length	Flow Rate
0.19 ml	10 cm	70 ml/min

A2010NWOR-PH

NeutrArt® Double Lumen Needlefree Valve Neutral, Luer Slip, Photoprotective

Priming Volume	Length	Flow Rate
0.42 ml	10 cm	70 ml/min

A3010NWOR-PH

NeutrArt® Triple Lumen Needlefree Valve Neutral, Luer Slip, Photoprotective

Priming Volume	Length	Flow Rate
0.57 ml	10 cm	70 ml/min