



Autofusion Enhance IV Sets




- *Autofusion Enhance* IV Sets with air vents and 15 μ m filters
- Hydrophilic membrane to prevent air entrainment if IV bag empties
- Available with 3-way connector and extension tubing
- Needle-free valves as optional extras


Autofusion Enhance IV Sets




1029DFAFE
Autofusion Enhance IV Set with air vent, 180cm tubing, rotating luer connector




1029DFVAFE
Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve




1029DFXAFE
Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector, 120cm extension tubing, rotating luer connector



1029DFXVAFE
Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, 120cm extension tubing, rotating luer connector



1029DFXMYAFE
Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector with needle-free valve, 120cm extension tubing, rotating luer connector



1029DFXVMYAFE
Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, with needle-free valve, 120cm extension tubing, rotating luer connector

Product code	Product description	Units per box	NHS Supply Chain code
1029DFAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, rotating luer connector	50	FSB2113
1029DFVAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve	50	FSB2114
1029DFXAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector, 120cm extension tubing, rotating luer connector	50	FSB2115
1029DFXVAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, 120cm extension tubing, rotating luer connector	50	FSB2116
1029DFXMYAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector, with needle-free valve, 120cm extension tubing, rotating luer connector	50	FSB2117
1029DFXVMYAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, with needle-free valve, 120cm extension tubing, rotating luer connector	50	FSB2450



Fannin

Caring for Life since 1829

www.fannin.eu

Fannin UK Ltd.
Fernhouse - Suite 3
Daniels Industrial Estate
Bath Road, Gloucestershire
GL5 3TJ, United Kingdom
T: +44 (0) 1453 874 130
E: salesll@fannin.eu

Fannin Ltd.
Fannin House
South County Business Park
Leopardstown
Dublin 18
T: +353 (0) 1290 7000
E: info@fannin.eu