

BARRIER surgical gown range folder



Whatever your priorities
BARRIER® surgical gowns have you covered

For the optimal balance in protection and comfort



Choosing the right surgical gown

is a balance between protection and comfort.

Surgical gowns are designed to protect patients and staff from cross-contamination, so protection is crucial when choosing a surgical gown. However, an uncomfortable surgical gown could put an entire procedure at risk – and comfort is a complicated thing.

Material that is breathable, soft, quiet and flexible enables staff to remain focused on the job at hand. And ergonomic elements improve the fit of the gown and provide greater freedom of movement.



Our sterile, procedure-specific products are folded for quick, easy and aseptic donning.

Mölnlycke Health Care has incorporated ergonomic elements into several* of our BARRIER® surgical gowns, including:

ERGONOMIC



Shoulder pleats



Forward facing sleeves




Long sleeves

* Ultimate, FPP and Classic

Choosing the right surgical gown for the right procedure.

Our range has been designed to provide the optimum combination of protection and comfort – from short, low-risk simple procedures to long, complicated high-fluid procedures. All BARRIER surgical gowns comply with both EN13795¹ and AAMI².




- PROTECTION REQUIREMENTS**
- Fluid levels
 - Duration

- COMFORT REQUIREMENTS**
- Breathability
 - Softness
 - Design and fit
 - A wide range of sizes


BARRIER surgical gown Ultimate

For the ultimate performance in protection and comfort

- Highly soft and noiseless material
- Maximum breathability for ultimate comfort
- Highly fluid repellent for advanced protection
- Ergonomic fit, for greater freedom of movement
- Available in Standard and High Performance gowns



ERGONOMIC




PROTECTION	COMFORT

BARRIER surgical gown Classic

For superior softness and comfort

- Superior softness for increased comfort
- Fluid repellent spunlace material for medium-risk protection
- Ergonomic fit, for greater freedom of movement
- Available in Standard and High Performance gowns



PROTECTION	COMFORT


ERGONOMIC



BARRIER surgical gown FFP*

For maximum fluid protection

- Lightweight, breathable material for thermal comfort
- Micropores combat perspiration, while polyethylene film prevents fluid penetration
- Impermeable to bodily fluids over the entire surface area
- Ergonomic fit, for greater freedom of movement



PROTECTION	COMFORT

ERGONOMIC




* Fluid Protection Plus, our most protective BARRIER surgical gown.

BARRIER surgical gown Universal

For all-round protection and simplicity

- Practical generous sizing for universal fit
- Ideal for interventions with medium to high amount of fluids and medium-duration procedures
- Long sleeves to prevent cuff movement
- Available in Standard and High Performance gowns




PROTECTION	COMFORT



BARRIER surgical gown Primary

For lightweight protection and simplicity

- Lightweight material for mobility and flexibility
- Fluid repellent for low-risk, short procedures
- T-cut design for good fit
- Available in Standard and High Performance gowns



PROTECTION	COMFORT



Surgical gowns from Mölnlycke Health Care

Mölnlycke Health Care has been leading the conversion to single-use surgical drapes and gowns for more than 35 years. As one of the most respected professional health care brands in the world, we are well known for our high quality products, customer-centric product development and innovative technologies.

Our unique range of products and services includes:

- Procedure-specific products for all surgical procedures.
- Consistently high levels of quality and safety, in our products and throughout our supply chain, always exceeding the strictest of standards.
- A commitment to reducing environmental impact and developing sustainable healthcare for the future.
- Lean supply and logistics solutions, in addition to expert training and support from our skilled and knowledgeable staff, to ensure you to make the most of your BARRIER surgical gowns.

At Mölnlycke Health Care, we understand the importance of finding the right balance between protection, comfort, efficiency and financial priorities. This is why we offer a wide range of carefully developed products, helping hospitals to achieve their ideal mix, maximise cost efficiency and free them up to treat more patients to better care.

1. European Standard EN 13795:2011 Surgical drapes, gowns and clean air suits, used as medical devices for patients, clinical staff and equipment – General requirements for manufacturers, processors and products, test methods, performance requirements and performance levels.

2. ANSI/AAMI PB70:2012 Liquid performance and classification of protective apparel and drapes for use in health care facilities.

3. FQN material study on competitor gown materials, SMS around 50-55 gsm, 2010.

Find out more at www.molnlycke.com

Mölnlycke Health Care AB, Box 13080, SE-402 52 Göteborg, Sweden.
T +46 31 722 30 00. F +46 31 722 34 00. www.molnlycke.com

The Mölnlycke Health Care and BARRIER® names and logos are trademarks registered globally to one or more of the Mölnlycke Health Care group of companies. Copyright (2014).