

BARRIER® Surgical mask

EXTRA PROTECTION

- 4-layer splash-resistant design
- >98% bacterial filtration efficiency
- Extra long tie bands
- EN14683:2005, type IIR
- Dispenser system supports efficient handling
- Visors protect from splashes
- Soft inner comfort layer
- Anti-fogging
- Cup Keeper reduces contact with mouth
- Ear-loops
- Tie-band

Proposed use: Wet surgery, when there is risk of splashes.



BARRIER® Surgical masks will help you treat more patients to better care

All BARRIER products from Mölnlycke Health Care are designed to offer the optimal combination of protection, functionality and comfort. Our innovative and procedure specific solutions come complete with logistic advantages and expert training and support. These elements together will create a ripple effect of savings in time, effort and costs - helping you treat more patients to better care.



658391



4224



4291



4232



4234



4236



4267



4313



4315



Dispenser system



Ref. no	Type	Tie-band	Ear-loop	Anti-fog	Visor	Splash resistant	Soft	Extra wide	Cup-keeper tie-band	Cup-keeper ear-loop	Anti-reflection tie-band	Anti-reflection ear-loop	Dispenser/Transport	Colour
658391	IIR	•				•	•		•				50/300	○
4224	IIR		•	•		•							50/300	●
4291	IIR	•		•		•							50/300	●
4232	IIR	•		•	•	•							25/100	●
4234	IIR	•		•		•							50/300	●
4236	IIR	•		•		•							50/300	●
4267	IIR		•	•	•	•							25/100	●
4313	IIR	•				•							50/600	●
4315	IIR		•			•							50/600	●