

BARRIER® Surgical mask

SPECIAL

- Low breathing resistance
- >98% bacterial filtration efficiency
- Extra long tie bands
- EN14683:2005, type II
- Dispenser system supports efficient handling

Proposed use: Dry surgery, when there is no risk of splashes.

- Visors protect from splashes
- Soft inner comfort layer
- Anti-fogging
- Cup Keeper reduces contact with mouth
- Anti-reflection
- Extra width for better facial coverage
- Ear-loops
- Tie-band





BARRIER® Surgical masks will help you treat more patients to better care

All BARRIER products from Mölnlycke Health Care are designed to offer the optimal combination of protection, functionality and comfort. Our innovative and procedure specific solutions come complete with logistic advantages and expert training and support. These elements together will create a ripple effect of savings in time, effort and costs - helping you treat more patients to better care.

BARRIER® Surgical mask

SPECIAL


















657170







658030

658200

658300 4229





4238





428801

Dispenser system

















		TIE-BAND	EAR-LOOP	ANTI-F06	VISOR	SPLASH RESISTANT	SOFT	EXTRA WIDE	CUP-KEEPER	CUP-KEEPER	ANTI-REFLECTION	ANTI-REFLECTION		
Ref. no	Туре	Tie-band	Ear-loop	Anti-fog	Visor	Splash resistant	Soft	Extra wide	Cupkeeper tie-band	Cupkeeper ear-loop	Anti Reflection ear- loop	Anti reflection tie-band	Dispenser/ Transport	Colour
657010		•		•									50/300	
657020		•		•								•	50/300	
657030		•		•	•							•	50/300	
657100		•						•					60/360	
657160		•						•					60/360	
657170		•						•					60/360	0
658000		•							•				60/360	
658010		•		•					•				50/300	
658020		•		•					•			•	50/300	
658030		•		•	•				•			•	50/300	
658200			•							•			50/300	
658300		•					•		•				50/300	0
4229			•	•									50/300	
4231		•		•									60/360	
4233		•											50/300	
4238		•		•									60/360	
428801		•		•									50/300	



31878

