Caressential

σcc Ý Ι Τ Α L

Caressential [™] Nitrile Examination Gloves are uniquely engineered to provide optimal performance and protection. Research was conducted based on user's needs to provide the most suitable product to cater to the requirements for protection, tactile sensitivity, feel and functionality. Users often require gloves with good tactile sensitivity to perform their task with confidence and providing the strongest film strength possible with ultimate tactile sensitivity.

Caressential [™] Examinations Gloves are certified with PPE Category III, providing the unsurpassed film strength and chemical resistance properties in the market. With proper contour design and surface treatment, the gloves are easy to don and provide the fit, feel and functionality desired by medical and speciality personnel.

Test Parameters

EN ISO 374: Determination of resistance to permeation by chemicals (EC Certification, PPE Category III) EN455, parts 1, 2 and 3 Medical Grade AQL 1.5 Examination Gloves Food Contact Test (21 CFR 177.2600) Chemotherapy Drugs Permeation Test (ASTM D 6978-05 (2013)) Skin Irritation & Sensitisation Test (ISO 10993-10:2010) Viral Penetration Test (ASTM F1671:13)

Code	Description	Case Size
D1700-22	Examination Gloves, Nitrile, Non-Sterile, Powder Free, Size: X-Small	10x200
D1701-22	Examination Gloves, Nitrile, Non-Sterile, Powder Free, Size: Small	10x200
D1702-22	Examination Gloves, Nitrile, Non-Sterile, Powder Free, Size: Medium	10x200
D1703-22	Examination Gloves, Nitrile, Non-Sterile, Powder Free, Size: Large	10x200
D1704-22	Examination Gloves, Nitrile, Non-Sterile, Powder Free, Size: X-Large	10x180





סכב עµד א דע דע דע דע דע דע דע

Nitrile Examination Gloves Freedom From Holes: AQL 1.5 Ambidextrous Beaded Cuff Finger Textured Powder-Free

To talk to our Fannin team about Caressential, please contact us at:



Fannin Ltd. Fannin House South County Business Park Leopardstown Dublin 18 T: +353 (0)1 290 7000 E: info@fannin.eu www.fannin.eu