



vernoncarus

First choice for **your healthcare partnership**

Detex™ X-ray detectable gauze swabs.

Detex™ X-ray detectable gauze swabs have provided Vernon-Carus with a presence in the operating room for over 4 decades, this market leading product has features that none of our smaller competitors have even considered.

Detex™ has a perforated inner pack and less packaging bulk than similar products. The comprehensive Detex™ swab range is manufactured to British Pharmacopia (BP) and EN 14079 standards ensuring a quality product in terms of strength, absorbency, thread construction and ability to be X-ray detectable.

Traceability is a crucial factor in the modern healthcare environment, all sterile Detex™ X-ray detectable gauze swabs have full traceability labels on all outer packs.

- Unique traceability system - provides 3 part recording in patient notes, operating room register and hospital records.
- Less Packaging - easier storage and helps toward a better environment.
- Perforated inner pack - makes presentation easier and faster.
- Automated production - provides consistence in quality and supply.
- Manufactured to BP and EN 14079 standards - ensures quality in absorbency, strength, construction and X- detectability.



detex®
X-RAY DETECTABLE GAUZE SWABS

Orders and Enquiries:

T +44 (0) 1772 299999 F +44 (0) 1772 299934
E vcorder@vernon-carus.co.uk

Vernon-Carus Ltd.

Matrix Park, Chorley, PR7 7NB, UK
T +44 (0)1772 299900 F +44 (0)1772 299901
E info@vernon-carus.co.uk W www.vernon-carus.co.uk



vernoncarus

First choice for **your healthcare partnership**

Detex™
X-ray detectable
gauze swabs.

Product description and size	Code	Case Quantity
5cm x 5cm 32ply	22482	96 pks
7.5cm x 7.5cm 8ply (Drug Tariff)	40603	144 pks
7.5cm x 7.5cm 8ply Green	22421	480 pks
7.5cm x 7.5cm 12ply	22414	480 pks
10cm x 7.5cm 32ply	28934	120 pks
10cm x 10cm 12ply Green	28931	160 pks
10cm x 10cm 16ply	28930	120 pks
10cm x 10cm 32ply	28933	80 pks
15cm x 2.5cm 16ply	21951	96 pks
15cm x 10cm 24ply	28932	72 pks
22.5cm x 22.5cm 12ply with tape	28937	36 pks
30cm x 30cm 12ply with tape	28938	24 pks
30cm x 30cm 12ply	28925	24 pks
30cm x 30cm 24ply with tape	28921	16 pks
36cm x 11cm 16ply	28936	32 pks
45cm x 45cm 6ply with tape	22410	48 pks
45cm x 45cm 12ply with tape	28939	12 pks

Orders and Enquiries:

T +44 (0) 1772 299999 **F** +44 (0) 1772 299934

E vcorder@vernon-carus.co.uk

Vernon-Carus Ltd.

Matrix Park, Chorley, PR7 7NB, UK

T +44 (0)1772 299900 **F** +44 (0)1772 299901

E info@vernon-carus.co.uk **W** www.vernon-carus.co.uk