



# Autofusion Enhance IV Sets




- *Autofusion Enhance* IV Sets with air vents and 15µm filters
- Hydrophilic membrane to prevent air entrainment if IV bag empties
- Available with 3-way connector and extension tubing
- Needle-free valves as optional extras

## Autofusion Enhance IV Sets




### 1029DFAFE

Autofusion Enhance IV Set with air vent, 180cm tubing, rotating luer connector




### 1029DFVAFE

Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve




### 1029DFXAFE

Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector, 120cm extension tubing, rotating luer connector




### 1029DFXVAFE

Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, 120cm extension tubing, rotating luer connector



### 1029DFXMYAFE

Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector with needle-free valve, 120cm extension tubing, rotating luer connector



### 1048B2MYAFE

Autofusion Enhance IV Set with air vent, 260cm tubing, back check valve, two needle-free valve Y-sites, rotating luer connector

Product code	Product description	Units per box
1029DFAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, rotating luer connector	50
1029DFVAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve	50
1029DFXAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector, 120cm extension tubing, rotating luer connector	50
1029DFXVAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, 120cm extension tubing, rotating luer connector	50
1029DFXMYAFE	Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector, with needle-free valve, 120cm extension tubing, rotating luer connector	50
1048B2MYAFE	Autofusion Enhance IV Set with air vent, 260cm tubing, back check valve, two needle-free valve Y-sites, rotating luer connector	50