



Fannin Longis Homecare Bed

Caring for Life since 1829



The **Fannin Longis Homecare Bed** combines safety and quality with comfort and good looks. Finished in attractive beech, the 4 Section Mattress Support and Bed Ends are equipped with individual Motor Systems. The Drive System complies with relevant standards and delivers uncompromising safety – allowing excellent head and foot sections adjustment. The Mattress Support Sections and Bed Ends are optimally designed for space-saving storage. Easily assembled by one person, the sturdy tubular frame provides stability and safety.

Features

- Solid Beech Head & Foot Boards
- Solid Beech Side Rails
- Separable for Transport
- 4 Section Mattress Support
- 2 Motor Adjustment
- Function Locking Splash Proof Handset
- Assistance Pole with Triangle
- Manual Position Adjustable Legrest
- Emergency Lowering Function
- Trendelenburg Option



Technical Data

■ Max Safe Working Load:	175kg
■ Max Patient Weight:	140kg
■ Total Product Weight:	102.5kg
■ Mattress Support Size:	90x200cm
■ Total Bed Dimensions:	102x215cm
■ Height Adjustment:	40cm to 80cm
■ Max Backrest Angle:	70°
■ Max Legrest Angle:	20°

Fannin Ltd.

Fannin House, South County Business Park
Leopardstown, Dublin 18, Ireland

T: +353 (0)1 290 7000 **E:** info@fannin.eu **W:** www.fannin.eu