

LateX Surgical Gloves Powder Free

Delivering Value:
Comfort, Precision and Reliable
Protection

FEATURES & BENEFITS

- **Quality Assured:** Every glove (100%) is air-inflation tested and visually inspected 7 times before dispatch
- **Skin-friendly formulation:** Our unique chlorine-free formulation production processes combined with the absence of chemical accelerators known to cause allergic responses mean skin irritations are minimised or eliminated. Tested according to ISO EN 10993
- **Easy donning:** Skin friendly polymer coated for easy donning even with damp hands saving you time.
- **Consistent calibrated grip:** The outer surface of Medispo surgical gloves are specially treated to provide a consistent grip. The grip is monitored and controlled during production and each lot must pass stringent internal grip controls before being released.
- **Complies with the following Standards and Regulations:**

Quality: ISO 9001, EN ISO 13485

Compliance EU: CE Mark, MDR, EN 455 1-4, EN ISO 374 1-5 Compliance USA: 510K, ASTM D3577

Environment: EN ISO 14001:2015

Product Information	
Latex / Latex-free/ material	Natural rubber latex
Cuff type	Beaded Cuff
Texture/ Grip	Micro Textured, Calibrated Grip, Medium
Powder / Powder-free	Powder-Free
Design	Anatomical, Forward Thumb Ball, Hand Specific
International Product Standards	EN 455 1-4, EN ISO 374 1-5:2016, EN ISO 21420:2019, ASTM D3577, ASTM D6978, ASTM F739,
Quality / Environmental Standards	ISO 9001, Iso 13485:2016, Iso 14001
Shelf Life	3 Years
Compliance to regulations	EU MDR; USA MDSAP
Notified body/ Test Laboratories	DEKRA, SATRA, ARDL, BSI
PPE compliance/ Type	PPE Regulation (EU) 2016/425, Type B
Biological Evaluation	EN 10993, Primary Skin Irritation Test, Dermal Sensitization Assay, Cytotoxicity Test
Additional tests	Viral Penetration PPE Compliance (ASTM F1671), Chemotherapy Drugs (ASTM D6978), Chemical Residue

Product Specification				
Item No.	Size	Length (±5mm)	Thickness (±0.03mm)	Pairs
MTPC-5.5	5.5#	280	Fingertip: 0.22	50/ Box
MTPC-6.0	6.0#	280		50/ Box
MTPC-6.5	6.5#	280	Palm: 0.19	50/ Box
MTPC-7.0	7.0#	285		50/ Box
MTPC-7.5	7.5#	285		50/ Box
MTPC-8.0	8.0#	290	Cuff: 0.16	50/ Box
MTPC-8.5	8.5#	290		50/ Box
MTPC-9.0	9.0#	300		50/ Box

* 4 dispenser boxes per shipping case, 200 pairs per shipping case

Find out more at www.fannin.eu

Fannin (UK) Ltd.

Westminster Industrial Est Repton Road, Measham Swadlincote,

Derbyshire DE12 7DT

T: +44 (0) 1530 514566

E: serviceuk@fannin.eu