



# OPERATING TABLE SHEET with bed transfer aid

- Skin-friendly, breathable, latex-free
- Maximum load: 200 kg
- Extremely tear resistant, lint-free, non-drip
- Comfosafe absorbent core reduces odours
- Proactive hygiene thanks to germ absorption
- Easy and safe repositioning
- Excellent price/performance ratio thanks to the low generation of waste
- **Transport and positioning with one pad**



# OPERATING TABLE SHEET with bed transfer aid



The operating table sheet with bed transfer aid is particularly absorbent, quickly absorbs fluids to protect the patient's skin and effectively prevents the risk of decubitus.

The heavy-duty, tear-resistant special carrier film of the patient pad can be used to take heavy patients weighing up to 150 kg safely to and from the operating table. The skin-friendly, smooth upper side is completely made of PP non-woven fabric. This ensures the surface remains dry. The hydrophobic, breathable lower side creates a germ barrier, providing additional hygiene that protects both the patients, yourself, your staff and your equipment from contamination.

Product Code	Outer dimensions	Core dimensions	Quantity absorbed (litres of fluid) max	Load	Packs
95929	210 x 80cm	200 x 60cm	2.5	200kg	1 PU = pack of 20