





















Suction products





SUCTION PRODUCTS


0G9019702	0G9029702	0G9039702	0G9049702	0G9059702
LINER	LINER	LINER	LINER	LINER
Receptal™ 2000mL Disposable Liner	Receptal™ 2000mL Disposable Liner with Valve	Receptal™ 1500mL Disposable Liner	Receptal™ 1500mL Disposable Liner with Valve	Receptal™ 2000mL Disposable Liner with Orthopaedic Port
				
				
Volume (L) 2	Volume (L) 2	Volume (L) 1.5	Volume (L) 1.5	Volume (L) 2
Nominal capacity (L) 1.9	Nominal capacity (L) 1.9	Nominal capacity (L) 1.4	Nominal capacity (L) 1.4	Nominal capacity (L) 1.9
Valved -	Valved Y	Valved -	Valved Y	Valved -
Tandem (non-valved) Y	Tandem (non-valved) -	Tandem (non-valved) Y	Tandem (non-valved) -	Tandem (non-valved) Y
Material PE/PP/PVC with DEHP	Material PE/PP/PVC with DEHP	Material PE/PP/PVC with DEHP	Material PE/PP/PVC with DEHP	Material PE/PP/PVC with DEHP
Customer Minimum Quantity Order (Eaches) 50	Customer Minimum Quantity Order (Eaches) 50	Customer Minimum Quantity Order (Eaches) 50	Customer Minimum Quantity Order (Eaches) 50	Customer Minimum Quantity Order (Eaches) 50
Receptal™ 2L Liner	Receptal™ 2L Liner	Receptal™ 1.5L Liner	Receptal™ 1.5L Liner	Receptal™ 2L Liner with 3/8inch patient port

0L6549702		0L6559702		0L6569702		0P5079702		0P5089702	
LINER		LINER		LINER		LINER		LINER	
Receptal™ 1500mL Disposable Liner with Pour Spout and Valve		Receptal™ 3000mL Disposable Liner with Valve. 3/8IN Port (0.95cm)		Receptal™ 3000mL Disposable Liner. 3/8IN Port (0.95cm)		Receptal™ 2000mL Disposable Liner with Pour Spout		Receptal™ 2000mL Disposable Liner with valve and Pour Spout	
									
									
Volume (L)	1.5	Volume (L)	3	Volume (L)	3	Volume (L)	2	Volume (L)	2
Nominal capacity (L)	1.4	Nominal capacity (L)	2.9	Nominal capacity (L)	2.9	Nominal capacity (L)	1.9	Nominal capacity (L)	1.9
Valved	Y	Valved	Y	Valved	-	Valved	-	Valved	Y
Tandem (non-valved)	-	Tandem (non-valved)	-	Tandem (non-valved)	Y	Tandem (non-valved)	Y	Tandem (non-valved)	-
Material	PE/ PP/ PVC with DEHP	Material	PE/ PP/ PVC with DEHP	Material	PE/ PP/ PVC with DEHP	Material	PE/ PP/ PVC with DEHP	Material	PE/ PP/ PVC with DEHP
Customer Minimum Quantity Order (Eaches)	50	Customer Minimum Quantity Order (Eaches)	50	Customer Minimum Quantity Order (Eaches)	50	Customer Minimum Quantity Order (Eaches)	50	Customer Minimum Quantity Order (Eaches)	50
Receptal™ 1.5L Liner with Pourspout		Receptal™ 3L Liner with 3/8inch patient port		Receptal™ 3L Liner with 3/8inch patient port		Receptal™ 2L Liner with Pourspout		Receptal™ 2L Liner with Pourspout	

100G9140610H3

ACCESSORIES

Suction Extension
Tubing Non-Sterile.
1.8m



Length (m)	1.8
------------	-----

Nominal capacity (L)	-
----------------------	---

Valved	-
--------	---

Tandem (non-valved)	-
---------------------	---

Material	PVC with DEHP
----------	---------------------


Customer Minimum Quantity Order (Eaches)	50
--	----

Suction tubing 1.8m
Receptal™

100L2060610H3

ACCESSORIES

Suction Tubing with
Finger. 1.8m



Length (m)	1.8
------------	-----

Nominal capacity (L)	-
----------------------	---

Valved	-
--------	---

Tandem (non-valved)	-
---------------------	---

Material	PVC with DEHP
----------	---------------------


Customer Minimum Quantity Order (Eaches)	50
--	----

1.8m non-sterile
suction tubing

100L3010540H3

ACCESSORIES

Receptal™ Gelling
Agent, 45g



Quantity (g)	45
--------------	----

Nominal capacity (L)	-
----------------------	---

Valved	-
--------	---

Tandem (non-valved)	-
---------------------	---

Material	PE / SP
----------	---------------


Customer Minimum Quantity Order (Eaches)	50
--	----

45g Bottle of Gelling Agent

100L860620763

ACCESSORIES

Gelling Agent, 30g



Quantity (g)	30
--------------	----

Nominal capacity (L)	-
----------------------	---

Valved	-
--------	---

Tandem (non-valved)	-
---------------------	---

Material	-
----------	---

Customer Minimum Quantity Order (Eaches)	50
--	----

30g Bag Gelling Agent

433036214

ACCESSORIES

Specimen Trap,
International. 4in
(10.4cm). Spout with
3/8in (0.95cm). Port,
Tethered Cap



Length (cm)	10.4
-------------	------

Nominal capacity (L)	-
----------------------	---

Valved	-
--------	---

Tandem (non-valved)	-
---------------------	---

Material	PE
----------	----

Customer Minimum Quantity Order (Eaches)	24
--	----

Specimen Trap Short

433036218

ACCESSORIES

Specimen Trap 8in (20.3cm). Spout with 1/2in. Port, Tethered Cap and Stiffener



Length (cm) 32

Nominal capacity (L) -

Valved -

Tandem (non-valved) -

Material PE

Customer Minimum Quantity Order (Eaches) 24

Specimen Trap Short

433039711

ACCESSORIES

Specimen Trap 4in. (10.15cm) Spout with 3/8in (0.95cm). Port, Tethered Cap and Stiffener



Length (cm) 32

Nominal capacity (L) -

Valved -

Tandem (non-valved) -

Material PE


Customer Minimum Quantity Order (Eaches) 24

Specimen Trap Short

433059220

ACCESSORIES

Attachment bracket for Receptal Accu-Measure Device for 1L, 1.5L and 2L Canisters



Volume (L) -

Nominal capacity (L) -

Valved -

Tandem (non-valved) -

Material PC


Customer Minimum Quantity Order (Eaches) 10

Plastic bracket for mounting Accumeasure™ to Canister

433066214

ACCESSORIES

Accu-Measure™ Receptal Accurate Measurement Device 250mL



Volume (L) -

Nominal capacity (L) -

Valved Y

Tandem (non-valved) -

Material PC


Customer Minimum Quantity Order (Eaches) 12

Accumeasure™

433096113

ACCESSORIES

VacGard™ External Automatic Shut-Off Valve



Volume (L) -

Nominal capacity (L) -

Valved Y

Tandem (non-valved) -

Material PC


Customer Minimum Quantity Order (Eaches) 10

VacGard™ external shut off valve

434559213

ACCESSORIES

Suction Extension
Tubing
Non-Sterile. 30m



Length (m) **30**

Nominal capacity (L) **-**

Valved **-**

Tandem (non-valved) **-**

Material **PVC
with
DEHP**


Customer Minimum Quantity Order (Eaches) **1**

30m PVC Suction Tubing

443009217

ACCESSORIES

Receptal™ Canister
with Tee Connector



Length (m) **-**

Nominal capacity (L) **-**

Valved **-**

Tandem (non-valved) **-**

Material **PC
with
Synthetic
Rubber
O
Ring**


Customer Minimum Quantity Order (Eaches) **5**

Receptal™ Canister T FITT Connector

M29205403

ACCESSORIES

Receptal™ Dump Box
2's



Quantity **2 liners**

Nominal capacity (L) **-**

Valved **-**

Tandem (non-valved) **-**

Material **Cardboard**


Customer Minimum Quantity Order (Eaches) **10**

Receptal™ Disposal Box for 2 used liners

M29305403

ACCESSORIES

Receptal™ Dump Box
6's



Quantity **6 liners**

Nominal capacity (L) **-**

Valved **-**

Tandem (non-valved) **-**

Material **Cardboard**


Customer Minimum Quantity Order (Eaches) **10**

Receptal™ Disposal Box for 6 used liners

434496114

CANISTER

Receptal™ Canister
1000 mL



Volume (L) **1**

Nominal capacity (L) **0.9**






Valved **-**

Tandem (non-valved) **-**

Material **Autoclavable
PC with
Synthetic
Rubber
O
Ring**

Customer Minimum Quantity Order (Eaches) **10**


Receptal™ 1L Canister

435036118		434456111		434238911		434259217		434329213	
CANISTER		CANISTER		CANISTER		HARDWARE		HARDWARE	
Receptal™ Canister 3000 mL		Receptal™ Canister 2000 mL		Receptal™ Canister 1500 mL		Receptal™ Canister Support, Short. 1000mL for Carousel or Wall Mounting (Use with Channeled Bracket)		Receptal™ Tandem Floor Stand with Conductive Castors, Stainless Steel	
									
Volume (L)	3	Volume (L)	2	Volume (L)	1.5	Volume (L)	-	Volume (L)	-
Nominal capacity (L)	2.9	Nominal capacity (L)	1.9	Nominal capacity (L)	1.4	Nominal capacity (L)	-	Nominal capacity (L)	-
Valved	-	Valved	-	Valved	-	Valved	-	Valved	-
Tandem (non-valved)	-	Tandem (non-valved)	-	Tandem (non-valved)	-	Tandem (non-valved)	-	Tandem (non-valved)	-
Material	Autoclavable PC with Synthetic Rubber O Ring	Material	Autoclavable PC with Synthetic Rubber O Ring	Material	Autoclavable PC with Synthetic Rubber O Ring	Material	Chrome plated	Material	Stainless Steel
Customer Minimum Quantity Order (Eaches)	8	Customer Minimum Quantity Order (Eaches)	10	Customer Minimum Quantity Order (Eaches)	10	Customer Minimum Quantity Order (Eaches)	5	Customer Minimum Quantity Order (Eaches)	1
Receptal™ 3L Canister		Receptal™ 2L Canister		Receptal™ 1.5L Canister		Receptal™ 1L Canister Support for wall and carousel mounting		Receptal™ 1.5-2L Tandem Canister Floor Stand on castor wheels	

434449217

HARDWARE

Channeled Wall Bracket for Wall Canister Support or Carousel Mounting



Volume (L) -

Nominal capacity (L) -

Valved -

Tandem (non-valved) -


Material **Stainless Steel**Customer Minimum Quantity Order (Eaches) **5**

Wall bracket for wall and carousel mounting

434479213

HARDWARE

Receptal™ in-line Quad Canister Stainless Steel Floor Stand, with Conductive Castors



Volume (L) -

Nominal capacity (L) -

Valved -

Tandem (non-valved) -


Material **Stainless Steel**Customer Minimum Quantity Order (Eaches) **1**

Receptal™ 1.5-2L Quad Canister Floor Stand on castor wheels

434269217

HARDWARE

Receptal™ Canister Support Long (Use with Channeled Bracket)



Volume (L) -

Nominal capacity (L) -

Valved -

Tandem (non-valved) -


Material **Chrome Plated**Customer Minimum Quantity Order (Eaches) **5**

Receptal™ 1.5-2L Canister Support for wall and carousel mounting

434299223

HARDWARE

Receptal™ Tandem Canister Support, Long (Use with Channeled Bracket)



Volume (L) -

Nominal capacity (L) -

Valved -

Tandem (non-valved) -


Material **Chrome Plated**Customer Minimum Quantity Order (Eaches) **1**

Receptal™ 1.5-2L Tandem Canister Support

434759212

HARDWARE

Tall Floor Stand with Small Head for Multiple Canisters, with Conductive Castors



Volume (L) -

Nominal capacity (L) -

Valved -

Tandem (non-valved) -

Material **Stainless Steel**Customer Minimum Quantity Order (Eaches) **1**

Receptal™ Carousel which holds a maximum of 4 Canisters