## ROBOOTER

REBOOT YOUR LIFE



# THE FIRST REBOOTER FOR BEGINEERS







## DBOOTER X40

#### STRUCTURAL ADVANTAGES

Well-designed, art-structured, ergonomic

- · 360° flexible large rocker
- · Independent shock suspension system
- · 3d ringback air seat
- Auti-decubitus memory foam cushion
- Rear anti-roll auxiliary wheel
- · Phone app remote control









Settings





Attitude Detection

Speed Regulation

### **SCENARIO SHOWS**

Simple operation, multi-scene use



**Turning easily** 

Omni-directional whell fit all kinds ofterrain

Putting into car trunk after folding

Removable lithium battery

Gull wing back light

PRODUCT PARAMETERS			
Endurance Mileage	ABOUT 25KM	Body Size	100.5 X 59 X93 CM
Max Speed	6KM/H	Fold Body Size	82 X 59 X 37 CM
Max Obstacle Height	6CM	Body Net Weight	ABOUT 41.8KG
Ditch Width	10CM	Charging Time	4-6H
Climbing Angle	10°	Front Tire	10-INCH OFF-ROAD OMNIDIRECTIONAL WHEEL
Chassis Net Height	7CM	Rear Tire	10 INCH PNEUMATIC GEAR TIRE
Minimum Rotation Radius	76CM	Side Arms	ANGLE ADJUSTABLE
Seat Height	57.5CM	Maximum Load	120KG
Seat Width	42CM	Battery	TERNARY LITHIUM BATTERY (20AH)
Seat Depth	42CM	Recommended Ambient Temperature	-10°C~50°C
Cushion	42 X 42CM	Waterproof Grade	IPX4