

The **Tulip GT** is a latex-free, multi-sized adult oropharyngeal airway designed to work in any environment



Tulip GT Oropharyngeal Airway

- Direct replacement for the Guedel and Facemask for all users
- 100% success first-contact
- Directly-connectable and Hands-free
- Increased O2 delivery >100% +
- Reduced regurgitation risk
- Self-extubation and In-line ET intubation



...it's easy to use for everybody!

Tulip GT Oropharyngeal Airway Features

No mask

- Reduces excess ventilating pressures and regurgitation risk
- Reduces leaks, increased O2 delivery

Hands free

- 100% success rates even with inexperienced users


Directly connectable

- 77%+ tidal volume and >100% + O2 delivery
- Ventilator connection
- Staff liberation – an airway for all users

Small size

- Reduced gagging for semi-conscious use
- Lowest cuff pressures known
- Low trauma

**Broad use for Scopes and Surgeries;
Bronchoscopy, Upper GI Endoscopy, ENT, Dental
Including Starved Elective surgery, NG Tube compatible**



Product Code	Description	Patient Size/ Sex	Patient approx. weight	Initial inflation volume & inflation pressure	Inflation pressure	Box Quantity
TGT-04	Tulip GT #4 wingless, red	Large Male	80-100kg approx.	+/-60mls	54.4cmH2O	25
TGT-03	Tulip GT #3 wingless, orange	Small Male/Large Female	60-80kg approx.	+/-50mls	54.4cmH2O	25
TGT-02	Tulip GT #2 wingless, green	Small Female/Large Child	40-80kg approx	+/-40mls	54.4cmH2O	25
TGT-01	Tulip GT #1 wingless, yellow	Large Child	40-60kg approx.	+/-40mls	54.4cmH2O	25
TGT-00	Tulip GT #0 wingless, blue	Medium Child	20-40kg approx.	+/-30mls	54.4cmH2O	25