



# **ANTIMICROBIAL Disposable Suction Liner Range**



## Our **commitment** is to **quality**









#### **Disposable Suction Liner Range**

The only system where all external surfaces are totally antimicrobial, designed to minimise the risk of cross infection.

Liners are available in 1 Litre, 2 Litre and 3 Litre versions. With canisters which fit any bracket within the hospital, suitable for ALL Clinical Areas.

Use as a completely closed system with the standard or PreGelled suction liners. Sample port versions for the addition of solidifying gel or for using with a sample sock are also available. All our liners with the exception of the cascade versions have a hydrophobic bacterial shut off filter for complete overflow protection.

Our unique patented designed liner is folded for ease of storage and begins to inflate immediately the vacuum is attached. With the integral carry handles, the liner allows for the safe removal from the canister, standing upright, whilst full of fluid, this makes

**BactiClear**<sup>®</sup> Antimicrobial Disposable Suction Range from VacSax<sup>®</sup> the product of choice





PART OF DCC

HARMA DEVICES LOGISTIC



#### VacSax Bacti**Clear®** Antimicrobial Suction Liner with Hydrophobic Bacterial Filters

Part No.	Description	Quantity
9910 <b>-</b> 340	1 Litre Antimicrobial Liner	box of 25
VAL-102	1 Litre Antimicrobial Wide Bore Liner	box of 30
VAL-201	2 Litre Antimicrobial Liner	box of 30
VAL-301	3 Litre Antimicrobial Liner	box of 20
VAL-202	2 Litre Antimicrobial Wide Bore Liner	box of 30
VAL-104	1 Litre Antimicrobial Liner with sample port	box of 30
VAL-204	2 Litre Antimicrobial Liner with sample port	box of 30
VAL-304	3 Litre Antimicrobial Liner with sample port	box of 20



### VacSax Bacti**Clear®** Antimicrobial Cascade Liners including Cascade/Tandem Tubing

Part No.	Description	Quantity
VAL-203	2 Litre Antimicrobial Cascade Liner including Cascade Tubing	box of 30
VAL-303	3 Litre Antimicrobial Cascade Liner including Cascade Tubing	box of 20
VAL-205	2 Litre Antimicrobial Sample Port Cascade Liner including Cascade Tubing	box of 30
VAL-306	3 Litre Antimicrobial Sample Port Cascade Liner including Cascade Tubing	box of 20







VacSax Bacti <b>Clear®</b> Antimicrobial PreGelled Liners with Hydrophobic Bacterial Filters				
Part No.	Description	Quantity		
9810 <b>-</b> 101	1 Litre PreGelled Antimicrobial Liner	box of 25		
GAL-211	2 Litre PreGelled Antimicrobial Liner	box of 30		
GAL-311	3 Litre PreGelled Antimicrobial Liner	box of 20		
GAL-212	2 Litre Wide Bore PreGelled Antimicrobial Liner	box of 30		



VacSax Bacti**Clear®** Antimicrobial PreGelled Cascade Liners including Cascade/Tandem Tubing

Part No.	Description	Quantity
GAL-213	2 Litre PreGelled Antimicrobial Cascad Liner including Cascade Tubing	e box of 30
GAL-313	3 Litre PreGelled Antimicrobial Cascad Liner including Cascade Tubing	e box of 20





Bacticlear Liners - Patent No. GB2497360

#### **VacSax Limited**

Western Wood Way, Langage Science Park, Plymouth, Devon PL7 5BG, U.K.

PART OF  $\mathbf{DCC}$  V

ITAL

PHARMA DEVICES LOGISTICS

Tel: +44 (0)1752 337000 Fax: +44 (0)1752 337600 Email: sales@vacsax.com

#### www.vacsax.com

