

# Fluent™

TRANSITIONLESS  
MICRO-INTRODUCER

**UP TO 57% LOWER**  
initial force than the leading competitor.\*



## Experience a **Fluent™** Introduction

Galt Medical Corp., one of the premier providers of .018" micro-access products, is pleased to introduce the Fluent™ Transitionless Micro-Introducer (TMI). The Fluent™ TMI is designed for ease of access, with redesigned angles and silicone coating. Test data shows the Fluent™ is more effective in reducing penetration forces than the leading competitor.

- **Fluent transition**
- **4Fr and 5Fr sizes**
- **Regular and Stiff cannula configurations**
- **X-Sharp® 21ga echogenic needle with bevel indicator**
- **Variety of guidewire options including high-radiopacity tips**
- **All-inclusive packaging**



\*1:1 simulated use test data on file at Galt Medical Corp.

For more information about the Fluent™ Transitionless Micro-Introducer or to set up a meeting with a Galt Sales Consultant call (800) 639-2800.

# Fluent™

## TRANSITIONLESS MICRO-INTRODUCER

### Fluent™ Transitionless Micro-Introducer Kit

Catalog Number	Sheath Size	Sheath Length	Guidewire	Guidewire Type	Tip Shape	Needle Length	Needle Type
KIT-038-01	4F	10cm	40cm	SS/Galt Gold®	Straight	7cm	Echo
KIT-038-02	5F	10cm	40cm	SS/Galt Gold®	Straight	7cm	Echo
KIT-038-03	4F	10cm	40cm	Nitinol/Galt Gold®	Straight	7cm	Echo
KIT-038-04	5F	10cm	40cm	Nitinol/Galt Gold®	Straight	7cm	Echo
KIT-038-05	4F	10cm	60cm	Nitinol/Galt Gold®	Straight	7cm	Echo
KIT-038-06	5F	10cm	60cm	Nitinol/Galt Gold®	Straight	7cm	Echo
KIT-038-07	4F	10cm	40cm	SS/SS	Straight	7cm	Echo
KIT-038-08	5F	10cm	40cm	SS/SS	Straight	7cm	Echo
KIT-038-17	4Fr	10cm	40cm	Nitinol/Galt Gold®	Straight	4cm	Echo
KIT-038-18	5Fr	10cm	40cm	Nitinol/Galt Gold®	Straight	4cm	Echo
KIT-038-19	4F	10cm	40cm	Nitinol/Galt Gold®, marked	Angled	7cm	Echo
KIT-038-20	5F	10cm	40cm	Nitinol/Galt Gold®, marked	Angled	7cm	Echo

### Stiffen Fluent™ Transitionless Micro-Introducer Kit

Catalog Number	Sheath Size	Sheath Length	Guidewire	Guidewire Type	Tip Shape	Needle Length	Needle Type
KIT-039-01	4F	10cm	40cm	SS/Galt Gold®	Straight	7cm	Echo
KIT-039-02	5F	10cm	40cm	SS/Galt Gold®	Straight	7cm	Echo
KIT-039-03	4F	10cm	40cm	Nitinol/Galt Gold®	Straight	7cm	Echo
KIT-039-04	5F	10cm	40cm	Nitinol/Galt Gold®	Straight	7cm	Echo
KIT-039-05	4F	10cm	40cm	SS/SS	Straight	7cm	Echo
KIT-039-06	5F	10cm	40cm	SS/SS	Straight	7cm	Echo
KIT-039-13	4F	10cm	60cm	Nitinol/Galt Gold®	Straight	7cm	Echo
KIT-039-14	5F	10cm	60cm	Nitinol/Galt Gold®	Straight	7cm	Echo
KIT-039-17	4Fr	10cm	40cm	Nitinol/Galt Gold®	Straight	4cm	Echo
KIT-039-18	5Fr	10cm	40cm	Nitinol/Galt Gold®	Straight	4cm	Echo
KIT-039-19	4F	10cm	40cm	Nitinol/Galt Gold®, marked	Angled	7cm	Echo
KIT-039-20	5F	10cm	40cm	Nitinol/Galt Gold®, marked	Angled	7cm	Echo

