

CliniMed Limited
 Freepost HY241
 High Wycombe
 Bucks
 HP10 8BR

2

To order Aura™ Drainable please call
 the CliniMed® Careline on 0800 036 0100

Aura™ Drainable pouch – order codes

Stoma size (mm)	Beige - Split cover			Clear		
	Mini	Midi	Maxi	Mini	Midi	Maxi
Cut-to-fit 13–60/80 (Mini 13–50/70)	MHDS 513	MHDM 513	MHDL 513	MHDS 713	MHDM 713	MHDL 713
Pre-cut 25	MHDS 525	MHDM 525	MHDL 525	-	-	-
Pre-cut 29	MHDS 529	MHDM 529	MHDL 529	-	-	-
Pre-cut 32	MHDS 532	MHDM 532	MHDL 532	-	-	-
Pre-cut 35	MHDS 535	MHDM 535	MHDL 535	-	-	-
Pre-cut 38	MHDS 538	MHDM 538	MHDL 538	-	-	-
Pre-cut 44	MHDS 544	MHDM 544	MHDL 544	-	-	-

Aura™ products are also available via NHS Supply Chain. Manufactured by Welland®, a CliniMed® Group company. Welland products are distributed in the UK by CliniMed Ltd. Tel: 01628 850100; Fax: 01628 523579. Email: enquiries@clinimed.co.uk or visit www.clinimed.co.uk. CliniMed Ltd, a company registered in England number 01646927. Registered office: Cavell House, Knaves Beech Way, Loudwater, High Wycombe, Bucks HP10 9QY. Welland®, CliniMed®, Aura™, HyperFlex®, Easy2peel™, Dual-Carb®, and SureFast™ are trademarks of CliniMed Holdings Ltd. ©2013 CliniMed Ltd. 1706/0213/1

CliniMed®



Manufactured by Welland®, a CliniMed® Group company



Drainable
Aura™
 with Manuka honey

page 8

CliniMed®

Our kindest drainable
 pouch yet



Drainable
Aura™
 with Manuka honey

page 1

NEW Aura™ Drainable

More freedom to live life with a stoma

New Aura™ Drainable features several innovations designed to give you the confidence and freedom to live life to the full.

- Modern, compact design in a range of sizes
- Tailored with pre-cut flat flange designs
- Fast, simple fitting and removal

New Aura™ Drainable pouches offer a comfortable, secure and discreet fit.



Manuka honey

One of the key new features of the Aura™ range is the introduction of medical grade **Manuka honey** to our skin-friendly Hyperflex® hydrocolloid. The inclusion of Manuka honey in our Hyperflex® hydrocolloid may help to promote healthy skin around the stoma.

Aura™ Drainable – our kindest pouch yet

- Unique HyperFlex® hydrocolloid flange with medical grade Manuka honey
- New Dual-Carb®i filter for better protection against leakage
- Easy-to-read flange cutting guide
- Secure integral outlet with large easy release tab
- Handy Easy2peel™ feature
- Super soft, water-repellent cover
- Choice of beige or clear, in three sizes
- Split back for easy application and monitoring of output



Request your free Aura™ sample

To receive a free sample of Aura™ Drainable, please fill in your details below and post the coupon

Name (Mr/Mrs/Miss/Ms): _____

Address: _____

Post Code: _____

Telephone No: _____ Mobile: _____

Email: _____

Details of the current pouch you use:

Name: _____ Code: _____ My exact stoma size is: _____ mm

Please send me:

Beige: Mini Midi Maxi Clear: Mini Midi Maxi

I understand that this request will be handled by CliniMed Limited or SecuriCare (Medical) Limited. I would like my details to be kept on file, so that I can be kept up to date with information about relevant new products and services.

1706/0213/1

page 2

page 7

Drainable
Aura™
with Manuka honey

Effective integral outlet


Trusted SureFast™ integral system

- SureFast™ integral system has an easy to open and reseal outlet for easy emptying and reduced cleaning time
- The large, easy release tab makes emptying, cleaning and closing simple and hygienic with no need for clips or ties
- Long outlet neck rolls up and is held securely with a strong hook and loop fastening system for a comfortable and secure leak-free fit
- Ideal for people with limited dexterity



“Best I’ve used in years” – Mrs MI, Derbyshire.

The benefits of Aura™ Drainable



page 6

page 3



Designed to give a comfortable and perfect fit

Fits snugly, protecting you

- Specifically designed to prevent leakage, the oval-shaped Hyperflex® hydrocolloid flange with Manuka honey is skin-friendly, satisfies adhesion and removal needs and is up to 22% larger than our previous flanges
- Extra adhesion at the 3 and 9 o'clock positions dramatically reduces separation and leakage in the area most at risk from pouch contents
- An improved cutting guide ensures a perfect fit around different stoma sizes and shapes, reducing skin irritation and potential leakage
- Curvilinear shape, contoured to match the shape of the body, ensures a snug, secure fit

Super soft and easy for you to use

- The super soft pouch cover offers exceptional comfort; and quick-drying fabric helps prevent perspiration build-up and makes life easier after bathing, showering or swimming
- A handy Easy2peel™ feature allows smooth removal of the pouch from the skin, making pouch changing quick and easy
- The split back allows easy application and monitoring, and a quiet film cover minimises rustling


"My skin looks very good now" – Mr EP, Dereham.

page 4

Advanced filter technology for better discretion

New improved Dual-Carb®i filter

- A new Dual-Carb®i filter, with a PTFE membrane, helps reduce leakage through the filter by creating a barrier against liquid
- Enhanced carbon formulation reduces odours, giving greater discretion over a longer time
- Laser cut filter allows maximum air flow management, helping to reduce ballooning
- The new Aura™ filter has advanced **deodorisation** and **permeability** when compared to a number of other pouches



Added filter protection

- Integrated film layer improves filter efficiency and reduces ballooning
- Seamless compact shape allows contents to move freely to bottom of pouch, reducing clearing time and pancaking

page 5

Drainable
Aura[™]
with Manuka honey