

# Single use *Integral* Silicone Laryngeal Masks

Standard & Reinforced ranges



- 100% medical-grade silicone
- Colour-coded inflation balloons for size identification
- Available in all sizes
- PVC-free
- Latex-free



Highly quality ISO connections



Highly flexible silicone tube



Silicone cuff enables atraumatic insertion & removal



Laryngeal mask introducer

### STANDARD *Integral Silicone Laryngeal Masks - single use*

Size	Product code	Units/box	Recommended patient weight	Inflation valve colour	NHS Supply Chain code
1	LM1004S2	10	Neonates/infants < 5 Kg	purple	FDD1000
1.5	LM1507S2	10	Infants 5-10 Kg	white	FDD1001
2	LM2010S2	10	Infants/children 10-20 Kg	red	FDD1002
2.5	LM2514S2	10	Children 20-30 Kg	orange	FDD1003
3	LM3020S2	10	Children 30-50 Kg	green	FDD1004
4	LM4030S2	10	Adults 50-70 Kg	yellow	FDD1005
5	LM5040S2	10	Adults 70-100 Kg	pink	FDD1006

### REINFORCED *Integral Silicone Laryngeal Masks - single use*

Size	Product code	Units/box	Recommended patient weight	Inflation valve colour	NHS Supply Chain code
1	LM1004S2R	10	Neonates/infants < 5 Kg	purple	FDD1007
1.5	LM1507S2R	10	Infants 5-10 Kg	white	FDD1008
2	LM2010S2R	10	Infants/children 10-20 Kg	red	FDD1009
2.5	LM2514S2R	10	Children 20-30 Kg	orange	FDD1010
3	LM3020S2R	10	Children 30-50 Kg	green	FDD1011
4	LM4030S2R	10	Adults 50-70 Kg	yellow	FDD1012
5	LM5040S2R	10	Adults 70-100 Kg	pink	FDD1013

### Laryngeal Mask Introducer - single use

Product code	Units/box	Recommended patient weight	NHS Supply Chain code
INT304	10	Medium 3 and 4	FDD1720

