




ARROW^gard BLUE[®] ACUTE HEMODIALYSIS CATHETER

In the ErgoPack™ System



RISK REDUCTION

INFECTION PROTECTION FEATURES:

- ARROW^gard Blue Acute Hemodialysis Catheter
- Procedural STOP Sign
- Central Line Insertion Checklist
- Central Line Insertion Practices Adherence Monitoring Form (CLIP Sheet)
- Hi-Lite Orange[®] ChloraPrep[®]
- Alcohol Hand Gel
- Sterile Hand Towel
- Maximal Barrier Set:
 - Surgical Cap
 - Mask with Face Shield
 - Surgical Gown
 - Exclusive Easy-to-use Head-to-toe Drape
- BioPatch[®]
- Tegaderm[®] Semi-permeable Transparent Dressing
- HemoHopper[®] Fluid Receptacle

SHARPS SAFETY FEATURES:

- SharpsAway™ II Locking Disposal Cup
- SafetyGlide[®] Injection Needles
- Safety Scalpel

The Acute Hemodialysis Catheter ErgoPack System provides 30 essential tools systematically organized for procedural ease. The ErgoPack System integrates the latest infection prevention and safety measures recommended by the CDC, SHEA/IDSA, IHI, INS, JCAHO, NHSN and OSHA into one convenient bundle.

To help reduce Central Line-Associated Bloodstream Infections (CLABSIs), the ARROW[®] Acute Hemodialysis Catheters utilize ARROW^gard, Teleflex's proprietary chlorhexadine-based antimicrobial technology. No other antimicrobial catheter technology has such a proven track record and broad-spectrum of antimicrobial protection. Since it was introduced in 1992, more than 30 studies and reviews have demonstrated that the ARROW^gard technology protects against the incidences of CLABSIs.¹

ARROW Acute Hemodialysis Catheters from Teleflex offer uncompromising quality to treat acute hemodialysis patients and safety features designed to protect against infection and the risk of vessel damage.






ARROW^gard Blue Acute Hemodialysis Catheter



The ErgoPack System

ERGOPACK SYSTEM ACUTE HEMODIALYSIS PRODUCTS*

CATHETER	PRODUCT CODE	DESCRIPTION	PER CASE
	CDC-22122-X1A	Arrow Large-Bore: 12 Fr. x 16 cm. 2 Lumens	5
	CDC-25122-X1A	Arrow Large-Bore: 12 Fr. x 20 cm. 2 Lumens	5
	CDC-26122-X1A	Arrow Large-Bore: 12 Fr. x 25 cm. 2 Lumens	5
	CDC-22142-X1A	Arrow Large-Bore: 14 Fr. x 16 cm. 2 Lumens	5
	CDC-25142-X1A	Arrow Large-Bore: 14 Fr. x 20 cm. 2 Lumens	5
	CDC-26142-X1A	Arrow Large-Bore: 14 Fr. x 25 cm. 2 Lumens	5
	CDC-23122-XU1A	Arrow You-Bend: 12 Fr. x 13 cm. 2 Lumens	5
	CDC-22122-XU1A	Arrow You-Bend: 12 Fr. x 16 cm. 2 Lumens	5
	CDC-25122-XU1A	Arrow You-Bend: 12 Fr. x 20 cm. 2 Lumens	5
	CDC-25142-XC1A	Arrow Precurved Large-Bore: 14 Fr. x 20 cm. 2 Lumens	5
	CDC-22142-XC1A	Arrow Precurved Large-Bore: Precurved: 14 Fr. x 16 cm. 2 Lumens	5

*For a customized option through our ARROW Select Kits program, contact your local Teleflex Dialysis Access Specialists.

REFERENCE:

1 <http://www.arrowintl.com/documents/pdf/literature/2008-0114.pdf>

Caution: US federal law restricts these devices to sale by or on the order of a physician.

TO PLACE AN ORDER CONTACT CUSTOMER SERVICE AT 1.800.523.8446

TELEFLEX MEDICAL Research Triangle Park, NC 27709
 Phone: 866.246.6990 Fax: 800.399.1028 Intl: 919.361.3961
TELEFLEX.COM

BioPatch is a registered trademark of Johnson and Johnson Corporation. Hi-Lite Orange and ChlorPrep are registered trademarks of CareFusion. SafetyGlide is a registered trademark of Becton, Dickinson, and Company. Tegaderm is a registered trademark of 3M Company. Arrow, ARROWgard Blue, ErgoPack, HemoHopper, SharpsAway, and Teleflex are trademarks or registered trademarks of Teleflex Incorporated or its affiliates. ©2011 Teleflex Incorporated. All rights reserved. 2011-0043 v1