



BARRIER® Surgical Gown - CLASSIC

- **Fluid repellent** spunlace for standard protection
- **Soft to the touch** and breathable for high comfort
- **Ergonomic design** for comfort and freedom of movement
- Suitable for all procedures
- **Treated to provide** increased repellency to low surface tension fluids.



BARRIER Surgical Gown - CLASSIC

	Ref. No.	Description	Size		Disp/Trp	
Standard Performance	STANDARD PERFORMANCE					
76	650101	SG SP, two towels, wrapped	М	121 cm	20/40	
() (at)	650102	SG SP, two towels, wrapped	L	121 cm	20/40	
	650104	SG SP, two towels, wrapped	L-L	131 cm	20/40	
	650108	SG SP, two towels, wrapped	XL	132 cm	18/36	
	650110	SG SP, two towels, wrapped	XL-L	+ 1 x L-L	18/36	
	650114	SG SP, two towels, wrapped	XXL-L	148 cm	18/36	
	650113	SG SP, two towels, wrapped	XXL-XL	161 cm	16/32	
High Performance	HIGH PERFORMANCE					
	660101	SG HP, two towels, wrapped	L	121 cm	18/36	
	660102	SG HP, two towels, wrapped	L-L	131 cm	18/36	
	660105	SG HP, two towels, wrapped	XL	132 cm	20/40	
	660106	SG HP, two towels, wrapped	XL-L	145 cm	16/32	
	660111	SG HP, two towels, wrapped	XXL-L	148 cm	16/32	
	660110	SG HP, two towels, wrapped	XXL-XL	148 cm	16/32	
	SG = Surgical Gown SP = Standard Performance HP = High Performance					
Impermeable front reinforcement to reinforcement						

through. • Breathable for

comfort.

- minimise risk of strike extended well above elbow to minimise risk of strike through.
 - Absorbing and soft for comfortable feeling.

Gown materials	Gown Body & Sleeve	Nonwoven pulp/polyester fibres			
	Front Reinforcement	Polyethylene plastic film (breathable)			
	Sleeve Reinforcement	Polyethylene plastic film, viscose/polyester nonwoven			
Product classification	EN 13795 Standard Performance and AAMI Level 2 & EN 13795 High Performance and AAMI Level 3				
Packaging	Vacuum single packed gowns provided in a double box packaging system, ensures perfect condition and				
	optimises storage space.				
Environment	At Mölnlycke Health Care we have a long-standing commitment to minimise the impact of our activities on the				
	environment, without compromising the efficiency and reliability of our products and services. Sound				
	environmental practice is a priority in our production and at our sites around the world.				

Find out more at www.molnlycke.co.uk

Mölnlycke Health Care Ltd, The Arenson Centre, Arenson Way, Dunstable, Bedfordshire, LU5 5UL. T 0870 6060766. www.molnlycke.co.uk

The Mölnlycke Health Care name and logo and BARRIER® and ProcedurePak® are registered trademarks of Mölnlycke Health Care AB. Copyright (2012)

