





# BARRIER Surgical Gown - CLASSIC

	Ref. No.	Description	Size	Disp/Trp
	<b>STANDARD PERFORMANCE</b>			
	650101	SG SP, two towels, wrapped	M 121 cm	20/40
	650102	SG SP, two towels, wrapped	L 121 cm	20/40
	650104	SG SP, two towels, wrapped	L-L 131 cm	20/40
	650108	SG SP, two towels, wrapped	XL 132 cm	18/36
	650110	SG SP, two towels, wrapped	XL-L + 1 x L-L	18/36
	650114	SG SP, two towels, wrapped	XXL-L 148 cm	18/36
	650113	SG SP, two towels, wrapped	XXL-XL 161 cm	16/32
	<b>HIGH PERFORMANCE</b>			
	660101	SG HP, two towels, wrapped	L 121 cm	18/36
	660102	SG HP, two towels, wrapped	L-L 131 cm	18/36
	660105	SG HP, two towels, wrapped	XL 132 cm	20/40
	660106	SG HP, two towels, wrapped	XL-L 145 cm	16/32
	660111	SG HP, two towels, wrapped	XXL-L 148 cm	16/32
	660110	SG HP, two towels, wrapped	XXL-XL 148 cm	16/32

**SG = Surgical Gown SP = Standard Performance HP = High Performance**

- Impermeable front reinforcement to minimise risk of strike through.
- Breathable for comfort.
- Impermeable sleeve reinforcement extended well above elbow to minimise risk of strike through.
- Absorbing and soft for comfortable feeling.

<b>Gown materials</b>	Gown Body & Sleeve Front Reinforcement Sleeve Reinforcement	Nonwoven pulp/polyester fibres Polyethylene plastic film (breathable) Polyethylene plastic film, viscose/polyester nonwoven
<b>Product classification</b>	EN 13795 Standard Performance and AAMI Level 2 & EN 13795 High Performance and AAMI Level 3	
<b>Packaging</b>	Vacuum single packed gowns provided in a double box packaging system, ensures perfect condition and optimises storage space.	
<b>Environment</b>	At Mölnlycke Health Care we have a long-standing commitment to minimise the impact of our activities on the environment, without compromising the efficiency and reliability of our products and services. Sound environmental practice is a priority in our production and at our sites around the world.	

Find out more at [www.molnlycke.co.uk](http://www.molnlycke.co.uk)

Mölnlycke Health Care Ltd, The Arenson Centre, Arenson Way, Dunstable, Bedfordshire, LU5 5UL.  
T 0870 6060766. [www.molnlycke.co.uk](http://www.molnlycke.co.uk)

The Mölnlycke Health Care name and logo and BARRIER® and ProcedurePak® are registered trademarks of Mölnlycke Health Care AB. Copyright (2012)