

Caressential



———— Clinical Disposables Range ————





CONTENTS

GOWNS	3
HEADWEAR	6
SUPPLEMENTARIES	7



The Fannin business is focussed on providing quality products and value added services that support the rapidly changing requirements of the healthcare sector.

Clinicians are looking for customised clinical solutions, cost savings and improved product availability. All of this is needed to deliver better services and quicker healing outcomes for their patients. Being one of the most established healthcare distributors in Ireland, Fannin is able to offer all of the above together with excellent services.

The Fannin range of non-sterile disposable products has been produced to the highest specifications, thereby ensuring that healthcare workers and patients are as comfortable as possible in every situation.

Working with Fannin, customers will receive quality products which are competitively priced and supported by an established company and team. Fannin is well placed to provide a strong partnership to support your organisation and its needs well into the future.

PATIENT GOWNS



CD00 Sleeveless Patient Gown

A sleeveless gown with an opening at the back, and neck/waist belts. Ideal for use as a patient examination gown and a visitor gown. CD00 is manufactured from 52g SMS nonwoven fabric.

Sizes: One Size

CD01 Long Sleeve Patient Gown

A long sleeve examination gown with an opening at the back, and neck/waist belts. Ideal for use as a patient examination gown, radiology gown and visitor gown, CD01 is manufactured from 52g SMS nonwoven fabric.

Sizes: One Size



CD02 Short Sleeve Patient Gown

A short sleeved gown with neck and waist belts, and an opening to the back. Made from 52g SMS nonwoven fabric, CD02 is ideal for use as an examination gown and a visitor gown.

Sizes: One Size

PATIENT GOWNS

CD02LM Full Wrap Patient Gown

An examination gown for use on patients required to change into a gown pre-procedure. This modesty gown is made from 52g SMS nonwoven fabric and can be worn to the front or to the back.

Sizes: One Size



CD17 Chest Wrap Patient Gown

A modesty gown for use on patients who are required to change into a gown pre-procedure. CD17 is made from 52g SMS nonwoven fabric, and can be worn to the front or to the back.

Sizes: One Size



GOWNS



CD19 Impervious Isolation Gown

An impervious isolation gown with cuffs for use in isolation, CSSD, theatre and for visitor use. Manufactured from PPSB 18 GSM coated PE13 GSM. Yellow in colour

Sizes: One Size

CD18 Thumb Up Gown

A thumb up gown to prevent cross contamination. Manufactured from LDPE film 27-28 GSM, this gown is ideal for use as a hospital gown, laboratory gown and as a visitor gown. Blue in colour.

Sizes: One Size



CD25 Isolation Gown

An isolation gown with latex-free elastic cuffs. This gown is made from PPSB 25g/m² and is Intended for use in isolation, CSSD, theatre and for visitor use. Blue in colour.

Sizes: One Size



HEADWEAR

CD4058 Surgical Cap

A surgeon's cap with tie, manufactured from spunbond, for use in theatre.

Sizes: One Size



CD1595 Face Mask

A 3ply facemask manufactured from Spunbond PPSB 14gsm and compliant with ASTM F2100-07 (equivalent 691014683 Type II).

Sizes: One Size



CD20 Face Shield

A face shield with a sponge foam head band and elastic strap. The shield is anti-fog with 0.10mm thickness. Latex free for use in theatre.

Sizes: One Size



SUPPLEMENTARY



CD07T Bed Sheet


A single-use, disposable bed sheet folded. For use on hospital beds. White in colour and manufactured from 30g Chemicalbond (Rayon).

Size: 140cm x 240cm

CD08 Clinical Sheet

A non-sterile, disposable clinical sheet for use as a utility sheet, tray inner and in CSSD. Manufactured from 30g Chemicalbond (Rayon) and available in two sizes.

Sizes: M (38cm x 41cm) L (30cm x 60cm)



CD03 Pillow Case

A non-sterile, disposable pillow case manufactured from 30g Chemicalbond (Rayon).

Size: 50cm x 76cm

SUPPLEMENTARY

CD09T Trolley Cover

A disposable sheet for covering examination trolleys. Manufactured from 30g Chemicalbond (Rayon) and white in colour.


Size: 80cm x 180cm



CD06 Cot Sheet

A large, general use sheet manufactured from nonwoven 35g chemicalbond (Rayon) in white with blue stripes.

Size: 55cm x 38cm



CD21 Disposable Curtain

A disposable curtain for use around hospital beds and in a ward setting. Manufactured from SMMS 100gsm fire retardant materials. Blue in colour.

Size: 7.5m width with 2m drop



SUPPLEMENTARY

CD04 Dry Patient Wipes

Economical dry patient wipes for every day patient hygiene and general use.
Manufactured from Chemicalbond 100% Rayon 35g/m2.

Size: Each wipe measures 20cm x 32cm



www.fannin.eu

Distribution in UK by:
Fannin (UK) Ltd,
42-46 Booth Drive
Park Farm South
Wellingborough
Northamptonshire
NN8 6GT

Fannin Ltd
Fannin House
South County Business Park
Leopardstown
Dublin 18
T: +353 (0) 1290 7000
E: info@fannin.ie

Product	Colour	Quantity	Code
Sleeveless Patient Gown	Blue	6 x 25 per box (150)	CD00
Long Sleeve Patient Gown	Blue	6 x 25 per box (150)	CD01
Short Sleeve Patient Gown	Blue	6 x 25 per box (150)	CD02
Full Wrap Patient Gown	Blue	4 x 25 per box (100)	CD02LM
Chest Wrap Patient Gown	Blue	4 x 25 per box (100)	CD17
Thumb Up Gown	Blue	75 per box (75)	CD18
Impervious Isolation Gown	Yellow	50 per box (50)	CD19
Isolation Gown	Blue	5 x 10 per box (50)	CD25
Surgical Cap	Blue	6 x 100 per box (600)	CD4058
Face Mask	Blue	6 x 50 per box (300)	CD1595
Face Shield	Clear	60 per box (60)	CD20
Bed Sheet	White	50 per box (50)	CD07T
Clinical Sheet	White	Medium: 100 x 40 per box Large: 100 x 20 per box	CD08M/L
Pillow Case	White	5 x 100 per box (500)	CD03
Trolley Cover	White	6 x 50 per box (300)	CD09T
Cot Sheet	Blue/White Striped	4 x 100 per box (400)	CD06
Disposable Curtain	Blue	10 per box (10)	CD21
Patient Dry Wipes	White	40 packs x 100 wipes per box	CD04
Endoscopy Pants	Blue	5 x 10 per box (50)	CD24M/L