




- Shoulder pleats allow freedom of movement for the upper body.
- Forward facing sleeves correspond with natural working position to allow comfort and freedom of movement.
- Long sleeve increases freedom of movement and reduces risk of glove sliding.

## BARRIER® Surgical Gown - FPP\*

- **Our highest protection gown**, impermeable to body fluids over the total surface area
- **Ultralight**, breathable material that draws moisture away from the body and allows perspiration vapour to escape, providing thermal comfort
- **Ergonomic design** for comfort and freedom of movement
- **Suitable** for procedures with high fluid levels.

\* FPP = Fluid Protection Plus

# BARRIER Surgical Gown - FPP

Fluid Protection Plus High Performance	Ref. No.	Description	Size		Disp/Trp
	670102	SG HP, two towels, wrapped	L	121 cm	32/64
	670103	SG HP, two towels, wrapped	L-L	131 cm	30/60
	670104	SG HP, two towels, wrapped	XL	132 cm	28/56
	670105	SG HP, two towels, wrapped	XL-L	145 cm	28/56
	670107	SG HP, two towels, wrapped	XXL-L	148 cm	26/52
	670106	SG HP, two towels, wrapped	XXL-XL	160 cm	28/56
	<b>SG = Surgical Gown HP = High Performance</b>				

<b>Gown materials</b>	Gown Body & Sleeve	Polyethylene film and polyester/polyethylene nonwoven
<b>Product classification</b>	EN 13795 High Performance and AAMI Level 3	
<b>Packaging</b>	Vacuum packed gowns provided in a double box packaging system, ensures perfect condition and optimises storage space.	
<b>Environment</b>	At Mölnlycke Health Care we have a long-standing commitment to minimise the impact of our activities on the environment, without compromising the efficiency and reliability of our products and services. Sound environmental practice is a priority in our production and at our sites around the world.	

\* FPP = Fluid Protection Plus

Find out more at [www.molnlycke.co.uk](http://www.molnlycke.co.uk)

Mölnlycke Health Care Ltd, The Arenson Centre, Arenson Way, Dunstable, Bedfordshire, LU5 5UL.  
T 0870 6060766. [www.molnlycke.co.uk](http://www.molnlycke.co.uk)

The Mölnlycke Health Care name and logo and BARRIER® and ProcedurePak® are registered trademarks of Mölnlycke Health Care AB. Copyright (2012)