

Standard E-Sac

The original Espiner Sac

DEPLOY

- using 5mm instruments

POSITION

- with our easy orientation tabs and black stitching to simplify identification of the Sac's mouth

RETRIEVE


- via laparoscopic incision reducing hospital stay, cost and complication




Allowing instrument
free extraction



One-piece tear
resistant design with
unique PUR coating
for maximum strength
and durability




Integrated
tail allowing
external control



Standard E-Sac Features:



- E** Comprehensive solution covering laparoscopic procedures from node biopsy to bowel resection.
- E** Minimally invasive reducing; length of stay, financial strain and improved patient experience.
- E** Our unique coated ripstop nylon is resistant to tearing if punctured.
- E** Includes X-ray opaque marker.
- E** Polyurethane coating impermeable to biological fluids.
- E** Elongated tail providing external control throughout procedure.
- E** Easily deploy using 5mm instrument.
- E** Instrument free extraction
– no claw forceps needed to complete.



New Product Code	Trocar (mm)	Diameter (mm)	Length (mm)	Volume (ml)	Number Per Box	GTIN-13 Code
EMP70-INS-10	10	45	100	60	10	5060358120404
EMP100-INS-10	10	55	160	200	10	5060358120428
EMP160-INS-5	10 & 12	102	250	900	5	5060358120442

If you would like to know more about our original & innovative range, please call us on **+44 (0)1275 878 801** or email **espiner@fannin.eu**
www.espinermedical.com