


# FX-PRO Mattresses

- Designed from the ground up for Community Care applications
- ➔ Improved specification for increased performance
- ➔ Ideal for Nursing Homes, Care Homes, Community Hospitals and Home Care use
- ➔ A range of professional products delivering ready made solutions from mattresses, underlays, overlays and cushions to cotside bumpers and safety mats
- ➔ Proven use in 1000's of Care Homes throughout the UK



- ➔ Standard mattress size:  
1930mm x 910mm x 150mm

# FX-PRO Mattress Features



# FX-PRO59 Mattress

- FX-PRO59 - COMFORT MATTRESS
- ➔ For general use to aid comfort
- ➔ Single block foam (SB33/120)
- ➔ Four-way turn design
- ➔ Comfort specification<sup>1</sup>
- ➔ Patient weight guide<sup>2</sup> - 185kg / 29 stone



Order Code: FX-PRO59

<sup>1</sup> Indication of risk category - not a prescriptive criteria

<sup>2</sup> Weight limits are to be used as a guide only and in conjunction with clinical judgements

# FX-PRO70 Mattress

- FX-PRO70 - 5" PROFILED MATTRESS
- ➔ Aids free flow of air and reduces surface tension
- ➔ Single side castellated foam (SB30/125)
- ➔ Two-way turn design
- ➔ Medium risk specification<sup>1</sup>
- ➔ Patient weight guide<sup>2</sup> - 159kg / 25 stone




Order Code: FX-PRO70

<sup>1</sup> Indication of risk category - not a prescriptive criteria

<sup>2</sup> Weight limits are to be used as a guide only and in conjunction with clinical judgements

# FX-PRO63 Mattress

- FX-PRO63 - PROFILED MATTRESS
- ➔ Aids free flow of air and reduces surface tension
- ➔ Single foam construction (SB33/120)
- ➔ Castellated foam both sides
- ➔ Four-way turn design
- ➔ High risk specification<sup>1</sup>
- ➔ Patient weight guide<sup>2</sup> - 216kg / 34 stone



Order Code: FX-PRO63


<sup>1</sup> Indication of risk category - not a prescriptive criteria

<sup>2</sup> Weight limits are to be used as a guide only and in conjunction with clinical judgements

# FX-PRO73 Mattress

Supporting

- FX-PRO73 GELTEX® inside MATTRESS
- ➔ Combination of CMHR and GELTEX® inside foams provides softness, support and breathability
- ➔ Layered construction with GELTEX® inside (SB40/170 + Geltex)
- ➔ Top layer of GELTEX® inside with supportive foam base
- ➔ Non turn (rotatable design)
- ➔ Very High risk specification<sup>1</sup>
- ➔ Patient weight guide<sup>2</sup> - 216kg / 34 stone



Order Code: FX-PRO73

<sup>1</sup> Indication of risk category - not a prescriptive criteria

<sup>2</sup> Weight limits are to be used as a guide only and in conjunction with clinical judgements

# GELTEX<sup>®</sup> inside

- WHY GELTEX<sup>®</sup> inside ?
- ➔ Incorporates the essential requirements into a single high performance material

BREATHABILITY

PRESSURE DISTRIBUTION

BODY SUPPORT

in a way you don't get with other materials,  
because of its unique formula with gel properties.

The breathable, open cell structure of the new GELTEX<sup>®</sup> inside technology, provides the patient with a unique comfortable feeling.

GELTEX<sup>®</sup> inside is not temperature sensitive and gently pushes back to provide the all-important patient support

Optimal body support is obtained with GELTEX<sup>®</sup> inside technology. Resilience is 5 times higher compared to visco-elastic foam.