



Reusable Optics, Disposable Blade.
Low cost, zero capital investment

Optics

- ✓ Reusable 50 intubations
- ✓ Fully isolated from the patient

Blades

- ✓ Disposable two sizes

Docking Station

- ✓ Recharges battery
- ✓ Displays service life



Reusable and disposable

- ✓ No cross contamination
- ✓ Low cost per use
- ✓ Zero capital investment
- ✓ Zero maintenance cost
- ✓ Can be used in MRI locations

Zero capital investment

- ✓ Full deployment
- ✓ Instant availability

SP

Product code	Description	ETT Size	Mouth Opening
A-011	Airtraq-SP Regular (size 3)	7.0 - 8.5	16mm
A-021	Airtraq-SP Small (size 2)	6.0 - 7.5	15mm
A-031	Airtraq-SP Paediatric (size 1)	4.0 - 5.5	11.5mm
A-041	Airtraq-SP Infant (size 0)	2.5 - 3.5	11mm
A-061	Airtraq-SP Adult Nasal	n/a	15mm
A-071	Airtraq-SP Double Lumen	35 - 41 Fr.	19mm

AVANT

Product code	Description	ETT Size	Mouth Opening
A-501	Airtraq Avant Optics		
A-511	Airtraq Avant Regular Blade (size 3)	7.0 - 8.5	17mm
A-521	Airtraq Avant Small Blade (size 2)	6.0 - 7.5	17mm
A-590	Airtraq Avant Docking Station		

CAMERA

Product code	Description
A-360	Airtraq Wi-Fi Camera

Where and when you need it



www.fannin.eu

Distribution in UK by:
Fannin (UK) Ltd
42 Kingfisher Court
Hambridge Road, Newbury
West Berkshire, RG14 5SJ, UK
T: +44 (0) 1635 550 100
E: serviceuk@fannin.eu

Fannin Ltd.
Fannin House
South County Business Park
Leopardstown
Dublin 18
T: +353 (0) 1290 7000
E: info@fannin.ie

FANUK14-3826



SP AVANT

Wi-Fi Camera



The most **simple** solution for Video Intubation



All in one
Ready to use
Fully disposable

Adults

- Adult regular 7.0 to 8.5 ETT
- Adult small 6.0 to 7.5 ETT
- Double lumen 35 to 41 Fr.
- Adult Nasal -



Children

- Paediatric 4.0 to 5.5 ETT
- Infant 2.5 to 3.5 ETT



Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc.

Wi-Fi Camera

Optics


- Video Recording
- 2 Axis Flip Screen
- Wi-Fi
 - PC
 - Phone
 - Tablet



Guided Video Intubation

Proven Technology. 1+ Million Intubations

- Improves C-L grade
- Facilitates ETT insertion
- Safe for the patient
- Easy to learn, simple to use



- No stylet required
- Allows intubation from any position
- Reduces C-spine movement
- Simplifies ETT exchange

New display options available soon

- Snap-on connector and App for Smartphones
- Integrated adapter to FOB cameras