



Single use Bougies

- Original and malleable Bougies with excellent memory
- Adult and paediatric sizes

Original Bougie



Malleable Bougie



Original Introducers

- Manufactured from low density polyethylene
- Sterile, single use
- Maintains shape for ease of insertion
- Calibrated so depth of insertion can be safely observed
- Coude Tip option
- Paediatric size has calibrated markings on both ends, enabling use as a Straight or Coude Tip Bougie
- Latex-free
- Non-PVC

Malleable Introducers

- Manufactured with an internal wire
- Sterile, single use
- Graduation marks every 5cm
- Latex free
- Non-PVC

Ordering details

Product code	Product description	Units per box	NHS Supply Chain code
Adult 15 Fr x 70 cm			
09-0212-70	Single use Adult Coude Tip Bougie	10	FTC443
09-0212-72	Single use Adult Straight Tip Bougie	10	FTC442
09-0213-92	Single use Adult Malleable Bougie	10	FTC446
Paediatric 10 Fr x 70cm			
09-0211-70	Single use Paediatric Coude/Straight Tip Bougie	10	FTC444
09-0213-82	Single use Paediatric Malleable Bougie	10	FTC447

