


CUFF-SAFE®

Endotracheal tubes

Cuff-Safe® endotracheal tubes supplied with AccuCuff™ Cuff Pressure Indicator



- Supplied with AccuCuff™ Cuff Pressure Indicator
- AccuCuff™ allows for continuous visual monitoring of cuff pressure, eliminating the need for regular spot-checking
- Negative pressure graduation on the AccuCuff™ indicates when the cuff is fully deflated, reducing the risk of trauma during removal
- Extensive range of standard, reinforced and preformed cuffed tubes available
- Manufactured from transparent PVC
- Tracheal Cuff Pressure Indicator 20-29cm H₂O
- High Volume, Low Pressure Cuff
- Clear connector – EGBAT compliant
- Murphy Eye
- Single Use
- Latex Free
- DEHP Free
- 5-year shelf life allows long term storage and reduces waste

Ordering details

Standard with AccuCuff™

Size	Code	Size	Code	Units per box
2.0mm	–			–
2.5mm	ETT-P32-25	6.5mm	ETT-P32-65	10
3.0mm	ETT-P32-30	7.0mm	ETT-P32-70	10
3.5mm	ETT-P32-35	7.5mm	ETT-P32-75	10
4.0mm	ETT-P32-40	8.0mm	ETT-P32-80	10
4.5mm	ETT-P32-45	8.5mm	ETT-P32-85	10
5.0mm	ETT-P32-50	9.0mm	ETT-P32-90	10
5.5mm	ETT-P32-55	9.5mm	ETT-P32-95	10
6.0mm	ETT-P32-60	10.0mm	ETT-P32-100	10

Reinforced with AccuCuff™

Size	Code	Size	Code	Units per box
2.0mm	–			–
2.5mm	–	6.5mm	ETT-J32-65	10
3.0mm	ETT-J32-30	7.0mm	ETT-J32-70	10
3.5mm	ETT-J32-35	7.5mm	ETT-J32-75	10
4.0mm	ETT-J32-40	8.0mm	ETT-J32-80	10
4.5mm	ETT-J32-45	8.5mm	ETT-J32-85	10
5.0mm	ETT-J32-50	9.0mm	ETT-J32-90	10
5.5mm	ETT-J32-55	9.5mm	ETT-J32-95	10
6.0mm	ETT-J32-60	10.0mm	ETT-J32-100	10

Oral, Preformed, South Facing with AccuCuff™

Size	Code	Size	Code	Units per box
2.0mm	–			–
2.5mm	–	6.5mm	ETT-K32-65	10
3.0mm	ETT-K32-30	7.0mm	ETT-K32-70	10
3.5mm	ETT-K32-35	7.5mm	ETT-K32-75	10
4.0mm	ETT-K32-40	8.0mm	ETT-K32-80	10
4.5mm	ETT-K32-45	8.5mm	ETT-K32-85	10
5.0mm	ETT-K32-50	9.0mm	ETT-K32-90	10
5.5mm	ETT-K32-55	9.5mm	ETT-K32-95	10
6.0mm	ETT-K32-60	10.0mm	ETT-K32-100	10

Nasal, Preformed, North Facing with AccuCuff™

Size	Code	Size	Code	Units per box
2.0mm	–			–
2.5mm	–	6.5mm	ETT-B32-65	10
3.0mm	ETT-B32-30	7.0mm	ETT-B32-70	10
3.5mm	ETT-B32-35	7.5mm	ETT-B32-75	10
4.0mm	ETT-B32-40	8.0mm	ETT-B32-80	10
4.5mm	ETT-B32-45	8.5mm	ETT-B32-85	10
5.0mm	ETT-B32-50	9.0mm	ETT-B32-90	10
5.5mm	ETT-B32-55	9.5mm	ETT-B32-95	10
6.0mm	ETT-B32-60	10.0mm	ETT-B32-100	10