



Double Chamber Blood Sets

- **DEHP-free**
- **200µm filtration**
- **4mm ID tubing**
- **20 drops/mL**
- **White flotation ball**
- **Available with 3-way connector and extension tubing**




Double Chamber Blood Sets (4mm ID)




7017DF4 NDC 212224

Double Chamber Blood Set with 200µm filter, 142cm tubing (4mm ID), rotating luer connector




7017DF4X

Double Chamber Blood Set with 200µm filter, 142cm tubing (4mm ID), 3-way connector, 142cm extension tubing, rotating luer connector



7017DF4XMY

Double Chamber Blood Set with 200µm filter, 142cm tubing (4mm ID), 3-way connector, needle free access port, 142cm extension tubing, rotating luer connector



7017DF4D

Dual Spike, Double Chamber Blood Set with 200µm filter, 142cm tubing (4mm ID), rotating luer connector

