

PAT-TUBE®

Endotracheal tube
with integrated
throat/pharyngeal pack

Pat-Tube® is a specialised endotracheal tube with an integrated throat/pharyngeal cuff that is used instead of a loose throat pack, eliminating the risk of accidental retention. It offers clinicians a safe, clinical alternative whenever a throat/pharyngeal pack is required



- High volume, low pressure throat/pharyngeal cuff eliminates the issue of throat/pharyngeal pack retention and reduces the risk of sore throat following surgery
- Inflatable throat/pharyngeal pack delivers a complete fluid seal, preventing surgical debris and blood entering the lower airway or oesophagus
- Supplied with two AccuCuff™ Cuff Pressure Indicators enabling continuous visual monitoring of both the pharyngeal pack and tracheal cuff pressure during use

- Tracheal Cuff Pressure Indicator 20-29cm H₂O
- Pharyngeal Pack Pressure Indicator 15-25cm H₂O
- Oral and nasal versions allow for a wide variety of clinical applications
- Paediatric sizes available in most variants
- Murphy Eye
- Single Use
- Latex Free
- DEHP Free



Tracheal Cuff Pressure Indicator 20-29cm H₂O



Pharyngeal Pack Pressure Indicator 15-25cm H₂O



Throat/Pharyngeal and Tracheal Cuffs

Ordering details

Standard Pat-Tube®

Size	Code	Units per box
4.5mm	30045*	10
5.0mm	30050*	10
5.5mm	30055*	10
6.0mm	60060	10
6.5mm	60065	10
7.0mm	60070	10
7.5mm	60075	10
8.0mm	60080	10
8.5mm	60085	10
9.0mm	60090	10

* Uncuffed tube

Reinforced Pat-Tube®

Size	Code	Units per box
4.5mm	40045*	10
5.0mm	40050*	10
5.5mm	40055*	10
6.0mm	10060	10
6.5mm	10065	10
7.0mm	10070	10
7.5mm	10075	10
8.0mm	10080	10
8.5mm	10085	10
9.0mm	10090	10

* Uncuffed tube

Oral, Preformed, South Facing Pat-Tube®

Size	Code	Units per box
6.0mm	50060	10
6.5mm	50065	10
7.0mm	50070	10
7.5mm	50075	10
8.0mm	50080	10
8.5mm	50085	10
9.0mm	50090	10

Nasal, Preformed, North Facing Pat-Tube®

Size	Code	Units per box
6.0mm	20060	10
6.5mm	20065	10
7.0mm	20070	10
7.5mm	20075	10