



Sterisets® Pre-filled Saline Syringes




- Easy to use
- Clear labelling
- Sterile packaging
- Cost effective
- Bar coded
- Anti reflux

Sterisets® Pre-filled Saline Syringes

Sterisets® is well known in the medical device industry. With the experience of manufacturing pre-filled syringes for many years, Sterisets has developed a new production line for the pre-filled syringes. The whole process is performed in a cleanroom ISO 5 environment and is fully automated. All the syringes are terminally sterilised for clinical usage.



The colour of the closure cap and label identify the volume of the saline more easily.




Colour coded label for easy recognition



Each inner box contains clear instructions for use

Customisation

Sterisets can produce bespoke sets that include the pre-filled saline syringes, for catheterisation or haemodialysis.



Ordering Information

| Product code | Product description | Closure cap | Sterility | Units per box |
|--------------|--|-------------|--------------------|---------------|
| 14365SNSF | Pre-filled Saline Syringe - 10mL with 3mL solution | ○ | Sterile content | 125 |
| 14366SNS | Pre-filled Saline Syringe - 10mL with 5mL solution | ● | Sterile content | 125 |
| 14363SNS | Pre-filled Saline Syringe - 10mL | ● | Sterile content | 125 |
| 14364SNS | Pre-filled Saline Syringe - 20mL | ○ | Sterile content | 60 |
| 14365SF | Pre-filled Saline Syringe - 10mL with 3mL solution | ○ | Externally sterile | 40 |
| 14366S | Pre-filled Saline Syringe - 10mL with 5mL solution | ● | Externally sterile | 40 |
| 14363S | Pre-filled Saline Syringe - 10mL | ● | Externally sterile | 40 |
| 14364S | Pre-filled Saline Syringe - 20mL | ○ | Externally sterile | 25 |