




# Total Intravenous Anaesthesia Sets






High flow distal Codan valve (up to 280 mL/min) and anti-siphon valves (to prevent free-flow and back-flow) for the syringe pumps included



Microbore tubing used in all Tivatek™ kits ensures a low priming volume and reduces drug wastage



Anti-kink bridge as additional safety feature



High specification non-kink tubing offers a reduced risk of occlusion – even with several tight knots in the infusion line, drugs continue to run at flow rates set



2-part TIVA Infusion Sets with integral needle-free access ports and back check valve now available

## Ordering information

### Total Intravenous Anaesthesia Sets - clinical dead space 0.68 mL

Product code	Product description	Priming volume (mL)	Units per box	NHS supply chain code
<b>2-way</b>				
TT202E	2-way 2.0m TIVA Infusion Set with high flow valve	2.2	50	
TT202FE	2-way 2.0m 2-part TIVA Infusion Set with 2 needle-free access ports and valve	2.2	50	
TT252FE	2-way 2.5m 2-part TIVA Infusion Set with 2 needle-free access ports and valve	3.0	50	FSB2374
TT302E	2-way 3.0m TIVA Infusion Set with high flow valve	3.0	50	
<b>3-way</b>				
TT223E	3-way 2.2m TIVA Infusion Set with high flow valve	4.3	50	
TT223ER	3-way 2.2m TIVA Infusion Set with high flow valve and red connectors	4.3	50	
TT253E	3-way 2.5m TIVA Infusion Set with high flow valve	4.9	50	
TT253FE	3-way 2.5m 2-part TIVA Infusion Set with 2 needle-free access ports and valve	5.7	50	FSB2375
TT303E	3-way 3.0m TIVA Infusion Set with high flow valve	5.7	50	
<b>4-way</b>				
TT224ER	4-way 2.2m TIVA Infusion Set with high flow valve and red connectors	6.4	25	
TT304E	4-way 3.0m TIVA Infusion Set with high flow valve	8.3	25	

## Manufactured in the UK



Caring for Life since 1829

[www.fannin.eu](http://www.fannin.eu)

**Fannin (UK) Ltd.**  
Westminster Industrial Estate  
Repton Road, Measham,  
Swadlincote, Derbyshire,  
DE12 7DT  
T: +44 (0) 1530 514566  
E: [serviceuk@fannin.eu](mailto:serviceuk@fannin.eu)

**Fannin Ltd.**  
Fannin House  
South County Business Park  
Leopardstown  
Dublin 18  
T: +353 (0) 1290 7000  
E: [info@fannin.eu](mailto:info@fannin.eu)