


THROMBOSIS 18 / POST OP 24

MEASURING POINTS



THROMBOSIS 18 / POST-OP 24 CIRCUMFERENCE MEASUREMENTS (cm)

	S		M		L		XL		XXL
	normal	max	normal	max	normal	max	normal	max	normal
g	45-54	53-62	53-62	61-70	61-70	69-78	69-78	77-86	77-86
c	30-36		35-41		40-46		45-51		50-56
b	19-21		22-24		25-27		28-30		31-33

COMPOSITION: Polyamide - Elastam

THROMBOSIS 18 LENGTH MEASUREMENTS (cm)

	SHORT	REGULAR	LONG
A-D	35-39	39-44	44-50
A-G	63-72	72-82	82-93

POST-OP 24 LENGTH MEASUREMENTS (cm)

	SHORT	LONG
A-D	34-39	39-45
A-G	62-72	72-84