



# Fannin Longis Low Bed

Caring for Life since 1829



The **Fannin Longis Low Bed** features a low level mattress platform of 23cm. The solid beech wood construction delivers comfort, safety and stability as well as a striking design. Offering the highest level of functionality, the 4-section mattress support and bed ends are equipped with individual motor systems. The drive system allows adjustment of the head and foot sections, and height adjustment of the bed ends – complying with the relevant standards delivering uncompromising safety. The **Fannin Longis Low Bed** can be transported and easily assembled by one person.

## Features

- Low Profiling Bed from 23cm to 63cm Height Adjustment
- Solid Beech Wood Construction
- 4 Section Mattress Support
- 2 Motor Adjustment (separable for transport)
- Function Locking Handset
- Manual Position Adjustable Legrest
- Removable Mattress Retainers
- Emergency Lowering Quick-action Bolt Mechanism
- Trendelenburg Option



## Technical Data

- |                          |              |
|--------------------------|--------------|
| ■ Height Adjustment:     | 23cm to 63cm |
| ■ Max Safe Working Load: | 175kg        |
| ■ Max Patient Weight:    | 140kg        |
| ■ Total Product Weight:  | 102.5kg      |
| ■ Mattress Support Size: | 90x200cm     |
| ■ Total Bed Dimensions:  | 102x215cm    |
| ■ Max Backrest Angle:    | 70°          |
| ■ Max Legrest Angle:     | 20°          |

### Fannin Ltd.

Fannin House, South County Business Park  
Leopardstown, Dublin 18, Ireland

**T:** +353 (0)1 290 7000 **E:** [info@fannin.eu](mailto:info@fannin.eu) **W:** [www.fannin.eu](http://www.fannin.eu)