

NEX D1 DRY[®] SURGICAL SCRUB



Ergonomic shape to allow a perfect grip
High flexibility to ensure accurate cleaning
Soft and kind to the skin
With nailcleaner for accurate cleaning under the nails
Advanced bristles system
Easy peel-off single wrapping
Single use

NEX
MEDICAL Surgical scrub

www.nexmedical.com

sales@nexmedical.it

Hand Hygiene is the basis of a correct and effective infection prevention in the hospital environment.

Surgical hand antiseptics, to destroy transient micro-organisms and inhibit the growth of resident micro-organisms, is routinely carried out before undertaking invasive procedures.

Antiseptics may reduce the risk of surgical site infections in patients.



NEX D1 DRY - sterile EO surgical brush/sponge

Description

NEX D1 DRY is a single use brush/sponge with nailcleaner
Without DETERGENT
Sterilized by ETHYLENE OXIDE (EO)
Single use
Expiration: 5 years

Indications

Surgical scrubbing of hands for antiseptic pre-op procedures.
NEX D1 DRY is compatible with any detergent or antiseptic cleansing solution.

Materials

Brush: PE medical grade
Sponge: polyurethane
Nailcleaner: PP medical grade
All components are LATEX free

Microbiological & skin irritation tests

Hemolysis test (ISO 10993-4:2008)
Cytotoxicity test (ISO 10993-5:2009)
Skin irritation test (ISO 10993-10:2002 / Amd 1:2006)

Manufacturer

NEX MEDICAL
20010 Casorezzo (MI) - ITALY
reg. PMC N. 26/2011/Off.299 / UNI EN ISO 9001:2008 / UNI EN ISO 13485:2004

Product Ref: NEX D1 DRY

Colour code: green

REF. D1 DRY

NEX D1 DRY[®] SURGICAL[®] SCRUB

STERILE EO

Spazzolino/spugna per pulizia mani con puliscunghie
Brush/sponge for handwashing, with nailcleaner
Brosse/éponge pour les soins des mains, avec cure-ongles
Cepillo/esponja para lavarse las manos, con limpiauñas
Espunja/escova, seca com limpa-unhas
Handwaschbürsten, mit Nagelreiniger
Børste/swamp, tør til håndvask med neglerenser
Handenborstel, droog met nagelreiniger

NEX
MEDICAL
www.nexmedical.it

Nex Medical, 20010 Casorezzo (MI) - ITALY