

NEX IODIO P2[®] SURGICAL[®] SCRUB



Ergonomic shape to allow a perfect grip
High flexibility to ensure accurate cleaning
Soft and kind to the skin
With nailcleaner for accurate cleaning under the nails
Advanced bristles system
Easy peel-off single wrapping
Single use

NEX
MEDICAL Surgical scrub

www.nexmedical.com

sales@nexmedical.it

Hand Hygiene is the basis of a correct and effective infection prevention in the hospital environment.

Surgical hand antisepsis, to destroy transient micro-organisms and inhibit the growth of resident micro-organisms, is routinely carried out before undertaking invasive procedures.

Antisepsis may reduce the risk of surgical site infections in patients.



NEX IODIO P2 PVP-I 7,5% surgical brush/sponge

Description

NEX IODIO P2 is a single use surgical hand scrub/brush with nailcleaner impregnated with Povidone Iodine antiseptic solution.

Indications

Antiseptic surgical scrub. Disinfection of the hands prior to surgery procedures or in the operating theatre.

Formulation Composition

Ingredients for 100g.:
Povidone Iodine g. 7,5
Coformulants and water q.s. to g. 100
free iodine 0,75%

Microbiological & skin irritation tests

Bactericidal:
EN 1040, pr-EN 12054 (Surgical handrub and handwash)
pr-EN 12791 (Surgical hand disinfection), EN 1650 (Candida Albicans)
Dermal irritation test and sensitization
(OECD guidelines 404, 410, ISO 10993)

Materials

Brush: PE medical grade Sponge: Polyurethane
Nailcleaner: PP medical grade
All components are LATEX free

Manufacturer

NEX MEDICAL
20010 Casorezzo (MI) - ITALY
reg. PMC N. 26/2011/Off.299 / UNI EN ISO 9001:2008 / UNI EN ISO 13485:2004

Product Ref: NEX IODIO P2
Colour code: brown

NEX IODIO P2[®] SURGICAL[®] SCRUB

Surgical brush/sponge with nailcleaner - Povidone iodine 7,5% (0,75% available iodine)
NEX IODIO P2 is a single use surgical brush/sponge impregnated with approximately 20ml of Povidone iodine antiseptic solution

Uses: surgical hand scrub, for cleansing hands and arms prior to surgery

Active ingredients: Povidone iodine 7,5% (equivalent to 0,75% available iodine)
Direct action towards bacteria gram+, gram- EN 12791, EN 12054, EN 1040

Directions: wet hands and forearms with water. Scrub for at least 3 minutes paying particular attention to the nails, cuticles and interdigital spaces. Rinse and repeat scrub. Discard brush-sponge. Rinse hands and forearms thoroughly and dry with sterile cloth.

Warnings For external use only. (R36) Irritating to eyes (S2) Keep out of reach of children (S25) Avoid contact with eyes and ears. (S26) In case of contact with eyes, rinse immediately with plenty of water and seek medical advice (S46) If swallowed, seek medical advice immediately and show this container or label

Ingredients 100 g contain : Povidone-Iodine g. 7,5 ; Coformulants and water q.s.to g. 100

Italian Ministry of Health registration PMC n°18975 (D.G.F.D.M. /405) Store at room temperature
Manufacturer's Address: Nex Medical s.a.s. - 20010 Casorezzo (MI) - ITALY Latex free
Made in Italy www.nexmedical.it info@nexmedical.it