



BARRIER[®] Surgical Gown - PRIMARY

- **Fluid repellent** SMS for adequate protection
- **Lightweight** for added comfort
- **T-cut** design for good fit.

BARRIER Surgical Gown - PRIMARY

	Ref. No.	Description	Size	Disp.
Standard Performance 	STANDARD PERFORMANCE			
	98000622	SG SP, two towels, wrapped	L 126 cm	42
	98000623	SG SP, two towels, wrapped	L-L 150 cm	42
	98000624	SG SP, two towels, wrapped	XL 126 cm	42
	98000625	SG SP, two towels, wrapped	XL-L 150 cm	42
High Performance 	HIGH PERFORMANCE			
	98000722	SG HP, two towels, wrapped	L 126 cm	30
	98000723	SG HP, two towels, wrapped	L-L 150 cm	28
	98000724	SG HP, two towels, wrapped	XL 126 cm	28
	98000725	SG HP, two towels, wrapped	XL-L 150 cm	26
SG = Surgical Gown SP = Standard Performance HP = High Performance				

- Front and sleeve reinforcement are impermeable to minimise risk of strike through
- Ultrasonic welded sleeves for increased repellency
- Comfort layer for added comfort.

Gown materials	Gown Body & Sleeve Front Reinforcement Sleeve Reinforcement.	Nonwoven polypropylene fibres Polyethylene plastic film, polyester nonwoven Polyethylene plastic film, polyester nonwoven.
Product classification	EN 13795 Standard Performance & EN 13795 High Performance.	
Packaging	Single packed gowns provided in a protective bag and transport box, ensures protection against dust, moisture and damage during transport, handling and storage.	
Environment	At Mölnlycke Health Care we have a long-standing commitment to minimise the impact of our activities on the environment, without compromising the efficiency and reliability of our products and services. Sound environmental practice is a priority in our production and at our sites around the world.	

Find out more at www.molnlycke.co.uk

Mölnlycke Health Care Ltd, The Arenson Centre, Arenson Way, Dunstable, Bedfordshire, LU5 5UL.
T 0870 6060766. www.molnlycke.co.uk

The Mölnlycke Health Care name and logo and BARRIER® and ProcedurePak® are registered trademarks of Mölnlycke Health Care AB. Copyright (2012)