




## 100% PLAIN VISCOSE SURGICAL STOCKINETTE

- Tubinette is a surgical stockinette used to hold dressings securely in place and protect clothing from ointments and creams



## HOW TUBINETTE WORKS

TubINETTE is used to hold dressings securely in place and to protect clothing from ointments and creams.

TubINETTE is made from 100% viscose and is intended for use where dressings have to be changed frequently.

## BENEFITS OF TUBINETTE

- Holds dressings comfortably and securely in place
- Quick and easy to apply, with or without applicator
- Range of eight sizes makes it easy to retain dressings anywhere on the body
- Particularly suitable for use with frequent dressing changes

## MATERIAL PROPERTIES

- 100% Viscose
- Natural rubber latex is not a constituent of TubINETTE or its packaging

## AREAS OF USE

TubINETTE holds dressings securely and comfortably in place, without constriction or compression.

The tubular construction makes application quick and easy, particularly when using a tubular bandage applicator.

The range of sizes available makes it possible to retain a dressing on any part of the body. The convenience and low cost of stockINETTES make them highly cost-effective.

### TUBINETTE ORDERING INFORMATION

Art. No	Size	Application	Pieces Per Case	Pip Code	NHS Code
2416	01	Fingers and toes	40	029-3654	EGJ333
2417	12	Bulky finger dressing	26	029-3662	EGJ335
2419	56	Adult limbs	30	029-3688	EGJ339
2420	78	Large adult limbs	30	029-3696	EGJ341
2421	T1	Small bodies, adult heads	12	029-3704	EGJ343
2422	T2	Large bodies	12	029-3712	EGJ345

Supplied individually cartoned in 20-metre rolls



Mölnlycke Health Care Ltd, Arenson Centre, Arenson Way, Dunstable, Bedfordshire LU5 5UL, UK. Tel. +44 [0] 870 6060766, [www.molnlycke.co.uk](http://www.molnlycke.co.uk)

The Mölnlycke Health Care name and logo, TubINETTE™ are registered trademarks of Mölnlycke Health Care AB.

©2012 Mölnlycke Health Care. All rights reserved.

UK134550611

