


CranioCurve™

Preformed Mesh



BIOMET®
MICROFIXATION
One Surgeon. One Patient.™

A titanium cranioplasty solution that's ahead of the curve.

Biomet Microfixation introduces CranioCurve Preformed Mesh as part of our comprehensive cranioplasty portfolio, for efficient coverage of cranial defects in multiple anatomic regions.



Features

- Offers 5 preformed mesh pieces for off-the-shelf convenience.
- 0.6mm grade II titanium for a balance of strength and flexibility.
- Region specific anatomic shape for efficient and aesthetic placement.
- Etched labels on the implant for easy recognition.
- Compatible with our 1.5 Neuro system for easy integration.

Preformed 0.6mm Mesh

01-9520	Occipital - OCC
01-9521	Left Parietal - LPS
01-9522	Right Parietal - RPS
01-9523	Left Parietal, Large- LPL
01-9524	Right Parietal, Large - RPL

CranioCurve Mesh Fit Options



To learn more about CranioCurve™ preformed mesh, please contact us at:
1520 Tradeport Drive • Jacksonville, Florida 32218-2480 • Tel 904.741.4400
Toll-Free 800.874.7711 • Fax 904.741.4500 • Order Fax 904.741.3059
www.biometmicrofixation.com