


Total IV Device Disinfection



SwabArt[®]

SwabArt[®] Disinfecting Cap for swabable needle-free valves is a single use disinfecting product that protects patients from infection. The product is impregnated with an absorbent sponge, saturated with 70% Isopropyl Alcohol (IPA).

- Provides active disinfection for up to 7 days if not removed
- SwabArt[®] is compatible with the EN ISO 80369-7-standard
- Compatible with all swabable needle-free valves
- All materials are sterile and latex-free
- Strip design allows hanging on drip poles
- Supplied in boxes of 25 strips of 10 disinfecting caps or as 250 individual disinfecting caps



“One CLABSI (Central Line Associated Bloodstream Infection) prevention method is disinfection of intravenous access points..... Use of a disinfectant cap on IV needleless connectors.....was associated with decreased CLABSI and costs.”

“Disinfectant cap use was associated with an estimated savings of almost \$300,000 per year in the hospital studied.”
 Merrill, Katreena Collette *et al.* *American Journal of Infection Control* 2014 **42** (12) 1274-1277

“...easy to use...”, “..With the alcohol-impregnated caps, it’s one less thing to think about.”, “Better Safety, Better Resource Utilization”
 Wawrzyniak, Marcelina M., *American Journal of Infection Control* 2014 **42** (6) S15-S16

Ordering details

Code	Product description	Form	Caps per box	NHSSC code
C100	SwabArt® Disinfecting Cap	Strips of 10 Caps	250	FSW1516
C200	SwabArt® Disinfecting Cap	Individual Caps	250	-



Fannin UK Ltd.
 Fernhouse - Suite 3
 Daniels Industrial Estate
 Bath Road, Gloucestershire
 GL5 3TJ, United Kingdom
 T: +44 (0) 1453 874 130
 E: salesll@fannin.eu

