



# *Autofusion Enhance IV Sets*




- *Autofusion Enhance* IV Sets with air vents and 15µm filters
- Hydrophilic membrane to prevent air entrainment if IV bag empties
- Available with 3-way connector and extension tubing
- Needle-free valves as optional extras


## Autofusion Enhance IV Sets




**1029DFAFE**  
Autofusion Enhance IV Set with air vent, 180cm tubing, rotating luer connector




**1029DF2AFE**  
Autofusion Enhance IV Set with air vent, 240cm tubing, rotating luer connector




**1029DFVAFE**  
Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve




**1029DFXAFE**  
Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector, 120cm extension tubing, rotating luer connector




**1029DFXVAFE**  
Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, 120cm extension tubing, rotating luer connector




**1029DFXMYAFE**  
Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector with needle-free valve, 120cm extension tubing, rotating luer connector



**1029DFXVMYAFE**  
Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, with needle-free valve, 120cm extension tubing, rotating luer connector



**1048BMYAFE**  
Autofusion Enhance IV Set with air vent, 210cm tubing, needle free access port, back check valve, rotating luer connector



**1048B2MYAFE**  
Autofusion Enhance IV Set with air vent, 260cm tubing, two needle free access ports, back check valve, rotating luer connector

| Product code  | Product description   | Units per box | NDC SKU |
|---------------|---|---------------|---------|
| 1029DFAFE     | Autofusion Enhance IV Set with air vent, 180cm tubing, rotating luer connector  | 50            | 230532  |
| 1029DF2AFE    | Autofusion Enhance IV Set with air vent, 240cm tubing, rotating luer connector  | 50            | -       |
| 1029DFVAFE    | Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve   | 50            | 244461  |
| 1029DFXAFE    | Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector, 120cm extension tubing, rotating luer connector   | 50            | -       |
| 1029DFXVAFE   | Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, 120cm extension tubing, rotating luer connector                         | 50            | 240890  |
| 1029DFXMYAFE  | Autofusion Enhance IV Set with air vent, 180cm tubing, 3-way connector, with needle-free valve, 120cm extension tubing, rotating luer connector                   | 50            | -       |
| 1029DFXVMYAFE | Autofusion Enhance IV Set with air vent, 180cm tubing, back check valve, 3-way connector, with needle-free valve, 120cm extension tubing, rotating luer connector | 50            | -       |
| 1048BMYAFE    | Autofusion Enhance IV Set with air vent, 210cm tubing, needle free access port, back check valve, rotating luer connector   | 50            | -       |
| 1048B2MYAFE   | Autofusion Enhance IV Set with air vent, 260cm tubing, two needle free access ports, back check valve, rotating luer connector                                    | 50            | -       |