






Introducer Bougies

Latex-free introducers are both firm and flexible for ease of insertion, calibrated for accurate entry depth and available in adult and paediatric sizes





Product code	Size	Product description	Units per box	NHS Supply Chain code
9-0211-70	10Fr x 70cm	Paediatric, Coudé Tip	10	FTC444
				
9-0212-70	15Fr x 70cm	Adult, Coudé Tip	10	FTC443
				
9-0212-72	15Fr x 70cm	Adult, Straight Tip	10	FTC442
				
9-0212-80	15Fr x 70cm	Adult, Coudé Tip (Ported)	10	FSM6153
				

- Low density polyethylene provides proper stiffness for ease of insertion
- Paediatric size features calibrated markings on both ends allowing it to be used as a straight or coudé tip
- Calibrated to ensure accurate distance insertion
- Paediatric size fits in 4mm to 6mm tubes
- Adult size fits in 6mm to 11mm tubes
- Adult ported variant available

Malleable Endotracheal Tube Introducers


- Combines benefits of an introducer and stylette into one product
- Paediatric size features calibrated markings on both ends, allowing it to be used as a straight or coude tip

Product code	Size	Product description	Units per box	NHS Supply Chain code
9-0213-82	10Fr x 70cm	Paediatric Malleable	10	FTC447
				
9-0213-92	15Fr x 70cm	Adult, Malleable	10	FTC446
				

Bougie-To-Go™ Endotracheal Introducers

- Innovative compact packaging designed specifically for use by pre-hospital professionals
- Low density polyethylene provides proper stiffness for ease of insertion
- Calibrated to ensure accurate distance insertion
- Adult size fits in 6mm to 11mm tubes


Unwraps without distorting introducer shape; keeps tip in correct insertion position



Designed specifically for use by pre-hospital professionals



Conveniently packaged to fit in paramedic bags, pockets and compact spaces

Product code	Size	Product description	Units per box	NHS Supply Chain code
9-0212-82	15Fr x 60cm	Adult, Coudé Tip	10	FSM6154
				



Caring for Life since 1829

www.fannin.eu

Fannin (UK) Ltd.
Westminster Industrial Estate,
Repton Road, Measham,
Swadlincote, Derbyshire,
DE12 7DT
T: +44 (0) 1530 514566
E: serviceuk@fannin.eu