



Baby Fisiojet

Aerosol Kit with pacifier and nebuliser



Infant Aerosol Kit

- Made of soft plastic materials
- Anatomical shape
- Wide dispensing hole with physiological angle for optimal nasal inhalation.
- Available in kit with nebuliser and connection tube
- Nebuliser tested with: Albuterol Sulphate, Ipratropium Bromide, Cromolyn Sodium
- Technical characteristics of the nebuliser: MMAD 1.9; GSD 2.8; breathable fraction <math>< 5\mu</math> 70%

Product code	Product description	Units per box
KM-119	Aerosol kit complete with tube, nebuliser and connecting tube	1



Soft Seal Masks

- Non-inflatable mask
- Pliable cushion



Product code	Product description	Units per box
VR1600	Neonate Soft Seal Mask	20
VR2600	Infant Soft Seal Mask	20
VR3600	Child Soft Seal Mask	20