

SunStim™ Pro

The SunStim™ Pro Peripheral Nerve Stimulator delivers reliable, adjustable current in one-second twitch pulses – it features a built-in train of four with four pulses per two seconds. This deluxe nerve stimulator, ergonomically designed with a digital screen for quick visual assessment of functionality, comes complete with a removable silicone case, ball electrodes, 9V battery and 6' lead wires with alligator clips, as well as an instruction manual.

SunStim™ Pro (8-1053-63)


- Double burst: Two 60-ms bursts of 50 Hz separated by 0.75s
- Twitch: 1 pulse per second
- Tetanus: 50 Hz (50 pulses/s) or 100 Hz (100 pulses/s)
- Built-in train of four: 4 pulses in a period of 2 seconds
- Stimulus pulse: square wave monophasic pulses (200- μ s duration)



Colourful Silicone Pro Wrap available in 9 colours to personalise and protect your device

Easy-to-read illuminated digital screen displays DBS, Twitch, Tetanus, TOF, current intensity and battery life

Clearly visible, finger-actuated front face membrane controls



Side switch to increase and decrease level of stimulation

Ergonomic, compact shape fits comfortably in hand

Removable silicone wrap enhances grip, helps protect device and comes in a variety of colours

Ordering information

Product code	Product description	Units per box
8-1053-63	SunStim™ Pro Peripheral Nerve Stimulator	1
8-1053-GR	Silicone Pro Wrap Grey	1
8-1053-BL	Silicone Pro Wrap Blue	1
8-1053-PL	Silicone Pro Wrap Purple	1
8-1053-FS	Silicone Pro Wrap Fuchsia	1
8-1053-PK	Silicone Pro Wrap Pink	1
8-1053-RD	Silicone Pro Wrap Red	1
8-1053-OR	Silicone Pro Wrap Orange	1
8-1053-GN	Silicone Pro Wrap Green	1
8-1053-TL	Silicone Pro Wrap Teal	1
8-1053-76	Ball Electrode and Hollow Pin	1
8-1051-78	6' Lead, Snap and Shielded Pin	1
8-1051-79	6' Lead, Clip and Shielded Pin	1
8-1053-87	Zippered Pouch 8.5" x 6"	1
8-1053-95	Pole Mount Bracket for 8-1053-63	1



8-1053-76

8-1051-78

8-1051-79

8-1053-95



Neoprene Zippered Pouch can hold and protect your nerve stimulator, manual, accessories and extra pro wraps

SunStim™ Plus

The SunStim™ Plus Peripheral Nerve Stimulator delivers reliable, adjustable current in one-second twitch pulses - it features a built-in train of four with four pulses per two seconds. This product comes complete with ball electrodes, 9V battery and 6' lead wires with alligator clips, as well as an instruction manual.



- Double burst: 2 60-ms bursts of 50 Hz separated by .75s
- Twitch: 1 pulse per second
- Tetanus: 50 Hz (50 pulses/s) or 100 Hz (100 pulses/s)
- Built-in train of four: 4 pulses in a period of 2 seconds
- Stimulus pulse: square wave monophasic pulses (200- μ s duration) period of 2 seconds
- Stimulator output current is fully adjustable from 0-70 mA
- Pole mount bracket enables easy storage and retrieval from IV pole
- Visual illumination for pulse activity and low battery life
- Lightweight and pocket-size for easy use and carry



8-1053-76



8-1051-78



8-1051-79



8-1053-85



Neoprene Zippered Pouch can hold and help protect your nerve stimulator, manual and accessories

Ordering information

Product code	Product description	Units per box
8-1053-62	SunStim™ Plus Peripheral Nerve Stimulator	1
8-1053-76	Ball Electrode and Hollow Pin	1
8-1051-78	Lead, Snap and Shielded Pin	1
8-1051-79	Lead, Clip and Shielded Pin	1
8-1053-85	Pole Mount Bracket for 8-1053-62	1
8-1053-87	Zippered Pouch 8.5" x 6"	1