



Peripheral Nerve Stimulator SunStim™ Plus



The Peripheral Nerve Stimulator SunStim™ Plus delivers reliable, adjustable current in one-second twitch pulses. Features include a built-in train of four with four pulses per two seconds. Ball electrodes, battery and lead wires with alligator clips all included.

Peripheral Nerve Stimulator SunStim™ Plus

- Stimulator output current is fully adjustable from 0-70 mA

- **SunStim™ Plus (8-1053-62) switch functions:**

- Double burst: 2 60-ms bursts of 50 Hz separated by 0.75s

- Twitch: 1 pulse per second

- Tetanus: 50 Hz (50 pulses/s) or 100 Hz (100 pulses/s)

- Built-in train of four: 4 pulses in a period of 2 seconds

- Stimulus pulse: square wave monophasic pulses (200-µs duration)



- LED (green) battery light indicates that the unit is on, flashes on low battery

- LED (red) flashes with audible indicator sounds each time a stimulus pulse is generated

- Lightweight: 0.25lbs including battery

- Standard battery: 9 volt alkaline



8-1053-76



8-1053-78, -79



8-1053-62

Peripheral Nerve Stimulator SunStim™ Plus

Product code	Product description	Units per box
8-1053-62	SunStim™ Plus Peripheral Nerve Stimulator	1

Accessories

Product code	Product description	Units per box
8-1053-76	Ball Electrode and Hollow Pin	1
8-1051-78	Lead 6' Snap and Shielded Pin	1
8-1051-79	Lead 6' Clip and Shielded Pin	1